



LAWRENCE
LIVERMORE
NATIONAL
LABORATORY

Sample Proficiency Test exercise

Armando Alcaraz, Hugh Gregg, Carolyn Koester

February 10, 2006

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

This work was performed under the auspices of the U.S. Department of Energy by University of California, Lawrence Livermore National Laboratory under Contract W-7405-Eng-48.

Sample Proficiency Test Exercise

The current format of the OPCW proficiency tests has multiple sets of 2 samples sent to an analysis laboratory. In each sample set, one is identified as a sample, the other as a blank. This method of conducting proficiency tests differs from how an OPCW designated laboratory would receive authentic samples (a set of three containers, each not identified, consisting of the authentic sample, a control sample, and a blank sample).

This exercise was designed to test the reporting if the proficiency tests were to be conducted. As such, this is not an official OPCW proficiency test, and the attached report is one method by which LLNL might report their analyses under a more realistic testing scheme. Therefore, the title on the report “Report of the Umpteenth Official OPCW Proficiency Test” is meaningless, and provides a bit of whimsy for the analyses and readers of the report.

H. Gregg
A. Alcaraz
C. Koester



**ORGANISATION FOR THE PROHIBITION
OF CHEMICAL WEAPONS**

**Report of the
Umpteenth Official OPCW
Proficiency Test**

Laboratory code: 05

Total number of pages: 44¹

¹ Total number of pages including cover page and all attachments

TABLE OF CONTENTS







	Page no.
Cover page	1
Table of contents	2
Summary	
Participating Laboratory	3
Quality Assurance/Quality Control (QA/QC)	5
Names and structures of all reported compounds	7
Analytical techniques: O1/05	8
O2/05	9
Sample preparation description/flowchart	10
Organic Sample O1	
Compound number 1	
Analytical technique 1: GC/MS-EI	12
Analytical technique 2: GC/MS-CI	15
Analytical technique 3: GC/dFPD	18
Compound number 2	
Analytical technique 1: GC/MS-EI	20
Analytical technique 2: GC/MS-CI	23
Analytical technique 3: GC/dFPD	26
Organic Sample O2	
Compound number 3	
Analytical technique 1: GC/MS-EI	28
Analytical technique 2: GC/MS-CI	31
Analytical technique 3: GC/dFPD	34
Compound number 4	
Analytical technique 1: GC/MS-EI	36
Analytical technique 2: GC/MS-CI	39
Analytical technique 3: GC/ dFPD	42
Comments	44

SUMMARY: PARTICIPATING LABORATORY

1. Participating laboratory

Laboratory code	05
Name of the laboratory/institute participating in the test	Lawrence Livermore National Laboratory
Contact person	Mr. Armando Alcaraz
Address	PO Box 808, M/S L-178 7000 East Avenue Livermore, CA 94551
Telephone number	925-423-6889
Fax number	925-423-9014
Email address	Alcaraz1@llnl.gov

2. Analysts and authentication

	Name	Title	Pages*	Date**	Signature**
1	Mr. Armando Alcaraz	Principal Investigator	All	Feb. 3, 2006	
2	Dr. Hugh Gregg	Co-PI, Senior Chemist	All	Feb. 3, 2006	
3	Dr. Carolyn Koester	Research Scientist	Prep & GC/MS	Feb. 3, 2006	
4	Dr. Phil Pagoria	Research Scientist	Synthesis	Feb. 3, 2006	
5	Dr. Robert Maxwell	Research Scientist	NMR analysis	Feb. 3, 2006	
6	Ms. Tuijauna Mitchell-Hall	QA Manager	QA/QC	Feb. 3, 2006	
7					
8					
9					

* Page numbers defining the responsibility area of the analyst;

** Date and signature of the responsible analyst;

SUMMARY: QUALITY ASSURANCE / QUALITY CONTROL (QA/QC)

1. Status of the laboratory (tick where applicable)

- ☒ Accreditation accepted:
Year 2001 Accreditation body: American Association for Laboratory Accreditation
Scope of accreditation: Chemical
- ☐ Accreditation planned/pending:
Target year _____ Accreditation body: _____
Scope of accreditation: _____
- ☐ Not accredited.

2. Quality system (tick where applicable)

- ☒ Described in a Quality Assurance Manual/Handbook. Quality system in accordance with:
☐ ISO 900_____, ☐ EN 4500_____, ☒ ISO Guide 17025, ☐ Other: _____
- ☐ No quality system. Please, fill in question number 3.

3. QA/QC Summary (Summary of the applied quality assurance and quality control (QA/QC) procedures concerning sample preparation, calibration, and analysis. Requested only from laboratories without a quality system).



THE AMERICAN
ASSOCIATION
FOR LABORATORY
ACCREDITATION

ACCREDITED LABORATORY

A2LA has accredited

LAWRENCE LIVERMORE NATIONAL LABORATORY

Livermore, CA


for technical competence in the field of

Chemical Testing

The accreditation covers the specific tests and types of tests listed on the agreed scope of accreditation. This laboratory meets the requirements of ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration Laboratories" and any additional program requirements in the identified field of testing.

Presented this 12th day of April 2004.




President
For the Accreditation Council
Certificate Number 1914-01
Valid to February 28, 2006

For tests or types of tests to which this accreditation applies,
please refer to the laboratory's Chemical Scope of Accreditation.

**American Association for Laboratory Accreditation**SCOPE OF ACCREDITATION TO ISO/IEC 17025-1999

LAWRENCE LIVERMORE NATIONAL LABORATORY
FORENSIC SCIENCE CENTER – OPCW PROJECT
7000 East Avenue Mailstop L-178
Livermore, CA 94550
Armando Alcaraz Phone: 925 423 6889

CHEMICAL

Valid To: February 28, 2006

Certificate Number: 1914-01

In recognition of the successful completion of the A2LA evaluation process, accreditation is granted to this laboratory to perform the following types of Qualitative Tests for Chemicals related to Chemical Warfare Convention (CWC) in unknown samples:

Sample Preparation BB-SP5, BB-SP6, BB-SP8, BB-SP9,
LL-SP1, LL-SP2, LL-SP3, LL-SP4

TestProceduresSpectroscopy

Nuclear magnetic resonance

BB-NMR1

Capillary zone electrophoresis / UV detection

LL-CE1

Gas chromatography / Fourier Transform
Infrared Spectrometry

BB-IR1

Chromatography

Gas chromatography / Element Specific Detectors

BB-GC1

Gas chromatography / Mass Spectrometry

BB-MS1

Liquid chromatography / Atmospheric pressure
Chemical Ionization / Mass Spectrometry

BB-MS4

Liquid Chromatography / Electrospray Ionization
Mass Spectrometry

LL-MS1

Chain of Custody for Laboratory

LL-CC1

Work Instructions for the Preparation of Test Samples
for OPCW Proficiency TestsQDOC/LAB/WI/ PT2
(Minus Section 11 Confirmation)





(A2LA Cert. No. 1914-01) Revised 04/15/2004

Page 1 of 1



SUMMARY: NAMES AND STRUCTURES OF ALL REPORTED CHEMICALS

Laboratory code: 05

Sam. code	Chem. no*	Chemical name	Chemical Abstract number	Chemical Structure	Molecular formula	Schedule number	Comments**
O1	1	Diethyl methylphosphonate	683-08-9		C ₅ H ₁₃ O ₃ P	2.B.04	
O1	2	Dicyclohexyl methylphosphonate	7040-53-1		C ₁₃ H ₂₅ O ₃ P	2.B.04	
O2	3	Diisopropyl methylphosphonate	1445-75-6		C ₇ H ₁₇ O ₃ P	2.B.04	
O2	4	Diisopropyl ethylphosphonate	1067-69-2		C ₈ H ₁₉ O ₃ P	2.B.04	
O3		No scheduled compounds					

* Chemical number defined by the participating laboratory and used throughout the report for the reported chemical.

** Explanation for the reporting of non-scheduled compounds, details can be added in the comment section of the report.

Note: There must be an unbroken chain of evidence linking each reported chemical to the original sample from which an aliquot was prepared and analyzed for the identification of this chemical.

SUMMARY: ANALYTICAL TECHNIQUES

Laboratory code: 05 Sample code(s): O1/05

Chemical number*	Chemical name	Chemical analysed as	Analytical technique	Method page no.	Aliquot name
1	Diethyl methylphosphonate	<input checked="" type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____	GC/MS-EI GC/MS-CI GC/dFPD	12 15 18	CW-1-159-1-O1 CW-1-159-1-O1 CW-1-159-1-O1
2	Dicyclohexyl methylphosphonate	<input checked="" type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____	GC/MS-EI GC/MS-CI GC/dFPD	20 23 26	CW-1-159-1-O1 CW-1-159-1-O1 CW-1-159-1-O1
		<input type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____			
		<input type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____			

* Chemical number defined by the participating laboratory (see Summary: Names and Structures of All Reported Chemical);

SUMMARY: ANALYTICAL TECHNIQUES

Laboratory code: 05 Sample code(s): O2/05

Chemical number*	Chemical name	Chemical analysed as	Analytical technique	Method page no.	Aliquot name
3	Diisopropyl methylphosphonate	<input checked="" type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____	GC/MS-EI GC/MS-CI GC/dFPD	28 31 34	CW-1-159-2-O2 CW-1-159-2-O2 CW-1-159-2-O2
4	Diisopropyl ethylphosphonate	<input checked="" type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____	GC/MS-EI GC/MS-CI GC/dFPD	36 39 42	CW-1-159-2-O2 CW-1-159-2-O2 CW-1-159-2-O2
		<input type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____			
		<input type="checkbox"/> original chemical <input type="checkbox"/> methylated <input type="checkbox"/> silylated <input type="checkbox"/> other: _____			


* Chemical number defined by the participating laboratory (see Summary: Names and Structures of All Reported Chemical);

SAMPLE PREPARATION DESCRIPTIONLaboratory code: 05 Sample code(s): 01/05, 02/05, 03/05**1. Sample preparation**

Sample/ Aliquot Code	Specification of Sample/ Type of Sample Preparation	Amount/ Volume	Sample Preparation Procedures	End Volume	Resulting Aliquot Code	Analytical Technique(s)
O1/05	Organic	0.2 mL	No sample workup.	200 µL	CW-1-159-1-O1	GC/MS-EI GC/MS-CI GC/dFPD
O2/05	None.	0.2 mL	No sample workup.	200 µL	CW-1-159-2-O2	GC/MS-EI GC/MS-CI GC/dFPD
O3/05	None.	0.2 mL	No sample workup.	200 µL	CW-1-159-3-O3	GC/MS-EI GC/MS-CI GC/dFPD

2. Additional information

--



Note: This flowchart is for visualization only; see the preceding sample preparation description page for sample aliquot numbers

GC-EI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 01/05 Chemical number: 1

Aliquot codes:

Sample: CW-1-159-1-01



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:		
Flow rate:	<input type="checkbox"/> ml/min	<input checked="" type="checkbox"/> 32 cm/s	
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow		
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Solvent delay time:	3 min	Scan range:	30-600 m/z
Electron energy:	70 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library spectrum:	Source : <input type="checkbox"/> OCAD (code:) <input type="checkbox"/> NIST <input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):
Comments:	

File : G:\DATA\TEST PT (UK) 2006\J0004290.D
Acquired 24 Jan 2006 18:10 Sample : CW-1-159-1-O1
Method : CW Misc info:



File : G:\DATA\TEST PT (UK) 2006\J0004293.D
Acquired: 25 Jan 2006 9:25 Sample : CW-CK-1-145-2
Method : CW Misc info: diethyl methylphosphonate




EI chromatograms supporting identification of compound **1**; TIC on left; EIC (m/z **125**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-1-O1** from **O1/05**, retention time **10.27** min.

Bottom: Chromatograms of authentic reference standard of **Diethyl methylphosphonate**, retention time **10.30** min.

File : D:\DATA\TEST PT (UK) 2006\J0004290.D
Acquired : 24 Jan 2006 18:10 using AcqMethod CW
Sample Name: CW-1-159-1-O1
Misc Info :



File : D:\DATA\TEST PT (UK) 2006\J0004293.D
Acquired : 25 Jan 2006 9:25 using AcqMethod CW
Sample Name: CW-CK-1-145-2
Misc Info : diethyl methylphosphonate



El mass spectra of:

Top: Compound **1** in Organic sample **O1/05**, aliquot **CW-1-159-1-O1**

Bottom: Authentic reference standard of **Diethyl methylphosphonate** corresponding to compound **1**
(MW: 152)

GC-CI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 01/05 Chemical number: 1

Aliquot codes:

Sample: CW-1-159-1-O1



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:		
Flow rate:	<input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s		
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow		
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Reaction gas:	<input type="checkbox"/> Methane <input type="checkbox"/> Isobutane <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other:		
Solvent delay time:	3 min	Scan range:	55-600 m/z
Electron energy:	235 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained <input type="checkbox"/> RT GC/MS-EI <input type="checkbox"/> RT GC/MS-CI
Comments:	

File : D:\DATA\TEST PT (UK) 2006\K0002913.D
Acquired: 24 Jan 2006 11:28 Sample : 1uL of CW-1-159-1-O1
Method : PCI-NH3CW.M Misc info: PCI-NH3, HP-5MS



File : D:\DATA\TEST PT (UK) 2006\K0002918.D
Acquired: 24 Jan 2006 18:36 Sample : 1uL of DEMP Std, 145-2
Method : PCI-NH3CW.M Misc info: PCI-NH3, HP-5MS




CI chromatograms supporting identification of compound **1**; TIC on left; EIC (m/z **153**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-1-O1** from **O1/05**, retention time **10.00** min.

Bottom: Chromatograms of authentic reference standard of **Diethyl methylphosphonate**, retention time **10.03** min.

File :D:\DATA\TEST PT (UK) 2006\K0002913.D
Acquired : 24 Jan 2006 11:28 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of CW-1-159-1-O1
Misc Info : PCI-NH3, HP-5MS



File :D:\DATA\TEST PT (UK) 2006\K0002918.D
Acquired : 24 Jan 2006 18:36 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of DEMP Std, 145-2
Misc Info : PCI-NH3, HP-5MS



CI mass spectrum of:

Top: Compound **1** in Organic sample **O1/05**, aliquot **CW-1-159-1-O1**

Bottom: Authentic reference standard of **Diethyl methylphosphonate** corresponding to compound **1**
(MW: **152**)

GAS CHROMATOGRAPHY TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 01/05 Chemical number: 1

Aliquot codes:

Sample: CW-1-159-1-01


ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890 GC dual FPD
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:
Flow rate:	<input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.
Injector temperature:	250 °C
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)
Detector:	<input type="checkbox"/> AED → Element(s) = <input type="checkbox"/> NPD <input checked="" type="checkbox"/> FPD → <input checked="" type="checkbox"/> P-mode <input type="checkbox"/> S-mode <input type="checkbox"/> Other:
Comments:	

IDENTIFICATION

Chemical identified as:	<input checked="" type="checkbox"/> Original Chemical <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source of Reference: <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library RI:	Measured RI of identified Chemical = RI OCAD = → OCAD Code =
Comments	


File : C:\DATA\TESTPT-1\FPD-0104.D
Acquired : 1/23/06 3:55:52 PM
Method name: CW.M
Sample Name: CW-1-159-1-01
Misc Info : OPCW sample



Compound: **1**
Sample: **O1/05**
Aliquot: **CW-1-159-1-O1**
Ret. time: **10.40 min.**

GC/dFPD (phosphorous trace only) chromatogram of Organic sample.

File : C:\DATA\TESTPT-1\FPD-0107.D
Acquired : 1/23/06 6:16:15 PM
Method name: CW.M
Sample Name: CW-CK-1-145-2
Misc Info : ~ 5 ng/uL diethyl methylphosphonate



Compound: **1**
Sample: **reference std.**
Ret. time: **10.40 min.**

GC/dFPD (phosphorous trace only) chromatogram of authentic reference standard **Diethyl methylphosphonate.**

GC-EI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): O1/05 Chemical number: 2

Aliquot codes:

Sample: CW-1-159-1-O1



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He	<input type="checkbox"/> N ₂	<input type="checkbox"/> H ₂ <input type="checkbox"/> Other:
Flow rate:	<input type="checkbox"/> ml/min	<input checked="" type="checkbox"/> 32 cm/s	
Flow control:	<input type="checkbox"/> Constant Pressure	<input checked="" type="checkbox"/> Constant Flow	
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Solvent delay time:	3 min	Scan range:	30-600 m/z
Electron energy:	70 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library spectrum:	Source : <input type="checkbox"/> OCAD (code:) <input type="checkbox"/> NIST <input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):
Comments:	

File : D:\DATA\TEST PT (UK) 2006\J0004290.D
Acquired 24 Jan 2006 18:10 Sample : CW-1-159-1-O1
Method : CW Misc info:



File : D:\DATA\TEST PT (UK) 2006\J0004292.D
Acquired: 25 Jan 2006 8:38 Sample : CW-CK-1-145-1
Method : CW Misc info: dicyclohexyl methylphosphonate




EI chromatograms supporting identification of compound **2**; TIC on left; EIC (m/z **97**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-1-O1** from **O1/05**, retention time **24.26** min.

Bottom: Chromatograms of authentic reference standard of **Dicyclohexyl methylphosphonate**, retention time **24.27** min.

File :D:\DATA\TEST PT (UK) 2006\J0004290.D
Acquired : 24 Jan 2006 18:10 using AcqMethod CW
Sample Name: CW-1-159-1-O1
Misc Info :



File :D:\DATA\TEST PT (UK) 2006\J0004292.D
Acquired : 25 Jan 2006 8:38 using AcqMethod CW
Sample Name: CW-CK-1-145-1
Misc Info : dicyclohexyl methylphosphonate



El mass spectra of:

Top: Compound **2** in Organic sample **O1/05**, aliquot **CW-1-159-1-O1**

Bottom: Authentic reference standard of **Dicyclohexyl methylphosphonate** corresponding to compound **2** (MW: **260**)

GC-CI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 01/05 Chemical number: 2

Aliquot codes:

Sample: CW-1-159-1-O1



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:		
Flow rate:	<input type="checkbox"/> ml/min	<input checked="" type="checkbox"/> 32 cm/s	
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow		
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Reaction gas:	<input type="checkbox"/> Methane <input type="checkbox"/> Isobutane <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other:		
Solvent delay time:	3 min	Scan range:	55-600 m/z
Electron energy:	235 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained <input type="checkbox"/> RT GC/MS-EI <input type="checkbox"/> RT GC/MS-CI
Comments:	

File : D:\DATA\TEST PT (UK) 2006\K0002913.D
Acquired 24 Jan 2006 11:28 Sample : 1uL of CW-1-159-1-O1
Method : PCI-NH3CW.M Misc info: PCI-NH3, HP-5MS



File : D:\DATA\TEST PT (UK) 2006\K0002919.D
Acquired: 24 Jan 2006 19:24 Sample : 1uL of DCMP Std, 145-1
Method : PCI-NH3CW.M Misc info: PCI, NH3, HP-5MS




CI chromatograms supporting identification of compound **2**; TIC on left; EIC (m/z **261**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-1-O1** from **O1/05**, retention time **23.87** min.

Bottom: Chromatograms of authentic reference standard of **Dicyclohexyl methylphosphonate**, retention time **23.88** min.

File :D:\DATA\TEST PT (UK) 2006\K0002913.D
Acquired : 24 Jan 2006 11:28 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of CW-1-159-1-O1
Misc Info : PCI-NH3, HP-5MS



File :D:\DATA\TEST PT (UK) 2006\K0002919.D
Acquired : 24 Jan 2006 19:24 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of DCMP Std, 145-1
Misc Info : PCI, NH3, HP-5MS



CI mass spectrum of:

Top: Compound **2** in Organic sample **O1/05**, aliquot **CW-1-159-1-O1**

Bottom: Authentic reference standard of **Dicyclohexyl methylphosphonate** corresponding to compound **2** (MW: **260**)

GAS CHROMATOGRAPHY TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 01/05 Chemical number: 2

Aliquot codes:

Sample: **CW-1-159-1-01**

ANALYSIS METHOD


Instrument Manufacturer and Type:	Agilent 6890 GC dual FPD
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:
Flow rate:	<input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.
Injector temperature:	250 °C
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)
Detector:	<input type="checkbox"/> AED → Element(s) = <input type="checkbox"/> NPD <input checked="" type="checkbox"/> FPD → <input checked="" type="checkbox"/> P-mode <input type="checkbox"/> S-mode <input type="checkbox"/> Other:
Comments:	

IDENTIFICATION

Chemical identified as:	<input checked="" type="checkbox"/> Original Chemical <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source of Reference: <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library RI:	Measured RI of identified Chemical = RI OCAD = → OCAD Code =
Comments	

File : C:\DATA\TESTPT-1\FPD-0104.D
Acquired : 1/23/06 3:55:52 PM
Method name: CW.M
Sample Name: CW-1-159-1-01
Misc Info : OPCW sample

->




Compound: **2**
Sample: **O1/05**
Aliquot: **CW-1-159-1-O1**
Ret. time: **24.09 min.**

GC/dFPD (phosphorous trace only) chromatogram of Organic sample.

File : C:\DATA\TESTPT-1\FPD-0106.D
Acquired : 1/23/06 5:29:25 PM
Method name: CW.M
Sample Name: CW-CK-1-145-1
Misc Info : ~5 ng/L dicyclohexyl methylphosphonate

->



Compound: **2**
Sample: **reference std.**
Ret. time: **24.10 min.**

GC/dFPD (phosphorous trace only) chromatogram of authentic reference standard **Dicyclohexyl methylphosphonate.**

GC-EI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 02/05 Chemical number: 3

Aliquot codes:

Sample: CW-1-159-2-O2



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:		
Flow rate:	<input type="checkbox"/> ml/min	<input checked="" type="checkbox"/> 32 cm/s	
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow		
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Solvent delay time:	3 min	Scan range:	30-600 m/z
Electron energy:	70 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library spectrum:	Source : <input type="checkbox"/> OCAD (code:) <input type="checkbox"/> NIST <input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):
Comments:	

File : D:\DATA\TEST PT (UK) 2006\J0004291.D
Acquired 24 Jan 2006 18:58 Sample : CW-1-159-2-O2
Method : CW Misc info:



File : D:\DATA\TEST PT (UK) 2006\J0004294.D
Acquired: 25 Jan 2006 10:12 Sample : CW-CK-1-145-3
Method : CW Misc info: diisopropyl methylphosphonate




EI chromatograms supporting identification of compound **3**; TIC on left; EIC (m/z **123**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-2-O2** from **O2/05**, retention time **11.45** min.

Bottom: Chromatograms of authentic reference standard of **Diisopropyl methylphosphonate**, retention time **11.48** min.

File :D:\DATA\TEST PT (UK) 2006\J0004291.D
Acquired : 24 Jan 2006 18:58 using AcqMethod CW
Sample Name: CW-1-159-2-O2
Misc Info :



File :D:\DATA\TEST PT (UK) 2006\J0004294.D
Acquired : 25 Jan 2006 10:12 using AcqMethod CW
Sample Name: CW-CK-1-145-3
Misc Info : diisopropyl methylphosphonate



El mass spectra of:

Top: Compound **3** in Organic sample **O2/05**, aliquot **CW-1-159-2-O2**

Bottom: Authentic reference standard of **Diisopropyl methylphosphonate** corresponding to compound **3** (MW: 180)

GC-CI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 02/05 Chemical number: 3

Aliquot codes:

Sample: CW-1-159-2-O2



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:		
Flow rate:	<input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s		
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow		
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Reaction gas:	<input type="checkbox"/> Methane <input type="checkbox"/> Isobutane <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other:		
Solvent delay time:	3 min	Scan range:	55-600 m/z
Electron energy:	235 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained <input type="checkbox"/> RT GC/MS-EI <input type="checkbox"/> RT GC/MS-CI
Comments:	

File : D:\DATA\TEST PT (UK) 2006\K0002914.D
Acquired 24 Jan 2006 12:40 Sample : 1uL of CW-1-159-2-O2
Method : PCI-NH3CW.M Misc info: PCI-NH3, HP-5MS



File : D:\DATA\TEST PT (UK) 2006\K0002917.D
Acquired: 24 Jan 2006 17:48 Sample : 1uL of DIMP Std, 145-3
Method : PCI-NH3CW.M Misc info: PCI-NH3, HP-5MS




CI chromatograms supporting identification of compound **3**; TIC on left; EIC (m/z **181**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-2-O2** from **O2/05**, retention time **11.16** min.

Bottom: Chromatograms of authentic reference standard of **Diisopropyl methylphosphonate**, retention time **11.19** min.

File :D:\DATA\TEST PT (UK) 2006\K0002914.D
Acquired : 24 Jan 2006 12:40 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of CW-1-159-2-O2
Misc Info : PCI-NH3, HP-5MS



File :D:\DATA\TEST PT (UK) 2006\K0002917.D
Acquired : 24 Jan 2006 17:48 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of DIMP Std, 145-3
Misc Info : PCI-NH3, HP-5MS



CI mass spectrum of:

Top: Compound **3** in Organic sample **O2/05**, aliquot **CW-1-159-2-O2**

Bottom: Authentic reference standard of **Diisopropyl methylphosphonate** corresponding to compound **3** (MW: 180)

GAS CHROMATOGRAPHY TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 02/05 Chemical number: 3

Aliquot codes:

Sample: CW-1-159-2-02


ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890 GC dual FPD
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:
Flow rate:	<input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.
Injector temperature:	250 °C
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)
Detector:	<input type="checkbox"/> AED → Element(s) = <input type="checkbox"/> NPD <input checked="" type="checkbox"/> FPD → <input checked="" type="checkbox"/> P-mode <input type="checkbox"/> S-mode <input type="checkbox"/> Other:
Comments:	

IDENTIFICATION

Chemical identified as:	<input checked="" type="checkbox"/> Original Chemical <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source of Reference: <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library RI:	Measured RI of identified Chemical = RI OCAD = → OCAD Code =
Comments	


File : C:\DATA\TESTPT-1\FPD-0105.D
Acquired : 1/23/06 4:42:39 PM
Method name: CW.M
Sample Name: CW-1-159-2-O2
Misc Info : OPCW sample



Compound: **3**
Sample: **O2/05**
Aliquot: **CW-1-159-2-O2**
Ret. time: **11.54 min.**

GC/dFPD (phosphorous trace only) chromatogram of Organic sample.

File : C:\DATA\TESTPT-1\FPD-0108.D
Acquired : 1/23/06 7:03:01 PM
Method name: CW.M
Sample Name: CW-CK-1-145-3
Misc Info : ~ 5 ng/uL diisopropyl methylphosphonate



Compound: **3**
Sample: **reference std.**
Ret. time: **11.56 min.**

GC/dFPD (phosphorous trace only) chromatogram of authentic reference standard **Diisopropyl methylphosphonate**.

GC-EI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): O2/05 Chemical number: 4

Aliquot codes:

Sample: CW-1-159-2-O2



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:		
Flow rate:	<input type="checkbox"/> ml/min	<input checked="" type="checkbox"/> 32 cm/s	
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow		
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Solvent delay time:	3 min	Scan range:	30-600 m/z
Electron energy:	70 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input checked="" type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library spectrum:	Source : <input type="checkbox"/> OCAD (code:) <input type="checkbox"/> NIST <input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):
Comments:	

File : D:\DATA\TEST PT (UK) 2006\J0004291.D
Acquired 24 Jan 2006 18:58 Sample : CW-1-159-2-O2
Method : CW Misc info:



File : D:\DATA\TEST PT (UK) 2006\J0004295.D
Acquired: 25 Jan 2006 10:59 Sample : CW-CK-1-146-3
Method : CW Misc info: diisopropyl ethylphosphonate




EI chromatograms supporting identification of compound **4**; TIC on left; EIC (m/z **111**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-2-O2** from **O2/05**, retention time **13.00** min.

Bottom: Chromatograms of authentic reference standard of **Diisopropyl ethylphosphonate**, retention time **13.01** min.

File :D:\DATA\TEST PT (UK) 2006\J0004291.D
Acquired : 24 Jan 2006 18:58 using AcqMethod CW
Sample Name: CW-1-159-2-O2
Misc Info :



File :D:\DATA\TEST PT (UK) 2006\J0004295.D
Acquired : 25 Jan 2006 10:59 using AcqMethod CW
Sample Name: CW-CK-1-146-3
Misc Info : diisopropyl ethylphosphonate



El mass spectra of:

Top: Compound **4** in Organic sample **O2/05**, aliquot **CW-1-159-2-O2**

Bottom: Authentic reference standard of **Diisopropyl ethylphosphonate** corresponding to compound **4**
(MW: **194**)

GC-CI-MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 02/05 Chemical number: 4

Aliquot codes:

Sample: CW-1-159-2-O2



ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890/5973 GC/MSD		
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:		
Flow rate:	<input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s		
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow		
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input type="checkbox"/> On Column <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.		
Injector temperature:	250 °C		
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane		
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm		
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)		
Reaction gas:	<input type="checkbox"/> Methane <input type="checkbox"/> Isobutane <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other:		
Solvent delay time:	3 min	Scan range:	55-600 m/z
Electron energy:	235 eV	Scan time:	0.7 s
Ionisation polarity:	<input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	Mass resolution:	0.7 u
Comments:			



IDENTIFICATION

Compound identified as:	<input checked="" type="checkbox"/> Original compound <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source : <input checked="" type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial
<input type="checkbox"/> Not compared to reference chemical or library spectrum:	Intense ions in spectrum are explained <input type="checkbox"/> RT GC/MS-EI <input type="checkbox"/> RT GC/MS-CI
Comments:	

File : D:\DATA\TEST PT (UK) 2006\K0002914.D
Acquired 24 Jan 2006 12:40 Sample : 1uL of CW-1-159-2-O2
Method : PCI-NH3CW.M Misc info: PCI-NH3, HP-5MS



File : D:\DATA\TEST PT (UK) 2006\K0002920.D
Acquired: 24 Jan 2006 20:12 Sample : 1uL of DIEP
Method : PCI-NH3CW.M Misc info: PCI, NH3, HP-5MS




CI chromatograms supporting identification of compound **4**; TIC on left; EIC (m/z **195**) on right.


Top: Chromatograms of Organic sample, aliquot **CW-1-159-2-O2** from **O2/05**, retention time **12.70** min.

Bottom: Chromatograms of authentic reference standard of **Diisopropyl ethylphosphonate**, retention time **12.71** min.

File :D:\DATA\TEST PT (UK) 2006\K0002914.D
Acquired : 24 Jan 2006 12:40 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of CW-1-159-2-O2
Misc Info : PCI-NH3, HP-5MS



File :D:\DATA\TEST PT (UK) 2006\K0002920.D
Acquired : 24 Jan 2006 20:12 using AcqMethod PCI-NH3CW.M
Sample Name: 1uL of DIEP
Misc Info : PCI, NH3, HP-5MS



CI mass spectrum of:

Top: Compound **4** in Organic sample **O2/05**, aliquot **CW-1-159-2-O2**

Bottom: Authentic reference standard of **Diisopropyl ethylphosphonate** corresponding to compound **4**
(MW: **194**)

GAS CHROMATOGRAPHY TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 05 Sample code(s): 02/05 Chemical number: 4

Aliquot codes:

Sample: CW-1-159-2-02


ANALYSIS METHOD

Instrument Manufacturer and Type:	Agilent 6890 GC dual FPD
Carrier gas:	<input checked="" type="checkbox"/> He <input type="checkbox"/> N ₂ <input type="checkbox"/> H ₂ <input type="checkbox"/> Other:
Flow rate:	<input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s
Flow control:	<input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow
Injection mode:	<input type="checkbox"/> Split → Split ratio = <input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.
Injector temperature:	250 °C
Column brand/phase:	Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane
Column Length x ID x Film thickness:	30 m x 0.25 mm x 0.25 µm
GC temperature programme:	40 °C (3 min), 8 °C/min, 300 °C (3 min)
Detector:	<input type="checkbox"/> AED → Element(s) = <input type="checkbox"/> NPD <input checked="" type="checkbox"/> FPD → <input checked="" type="checkbox"/> P-mode <input type="checkbox"/> S-mode <input type="checkbox"/> Other:
Comments:	

IDENTIFICATION

Chemical identified as:	<input checked="" type="checkbox"/> Original Chemical <input type="checkbox"/> Methyl ester derivative <input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative <input type="checkbox"/> TMS (Trimethylsilyl) derivative <input type="checkbox"/> Other derivative:
Retention parameter used for (peak) identification:	<input checked="" type="checkbox"/> Retention time (Rt) <input type="checkbox"/> Scan number
<input checked="" type="checkbox"/> Compared to reference chemical:	Source of Reference: <input checked="" type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial
<input type="checkbox"/> Compared to library RI:	Measured RI of identified Chemical = RI OCAD = → OCAD Code =
Comments	


File : C:\DATA\TESTPT-1\FPD-0105.D
Acquired : 1/23/06 4:42:39 PM
Method name: CW.M
Sample Name: CW-1-159-2-O2
Misc Info : OPCW sample



Compound: **4**
Sample: **O2/05**
Aliquot: **CW-1-159-2-O2**
Ret. time: **13.04 min.**

GC/dFPD (phosphorous trace only) chromatogram of Organic sample.

File : C:\DATA\TESTPT-1\FPD-0109.D
Acquired : 1/24/06 3:15:33 PM
Method name: CW.M
Sample Name: CW-CK-1-146-3
Misc Info : diisopropyl ethylphosphonate



Compound: **4**
Sample: **reference std.**
Ret. time: **13.04 min.**

GC/dFPD (phosphorous trace only) chromatogram of authentic reference standard **Diisopropyl ethylphosphonate**.

COMMENTS

1. General

No relevant CWC schedule compounds were identified in sample O3/05

2. Sample preparation

None

3. Analysis

No additional comments

4. Report

The field length in the Chemstation software often did not allow for the full chemical name and may be truncated or abbreviated.

The GC Chemstation report software often truncated the time axis at approximately 80% of full scale; in no case was data of significance not displayed correctly.