



LAWRENCE  
LIVERMORE  
NATIONAL  
LABORATORY

# Pentaerythritol Tetranitrate (PETN) Surveillance by HPLC-MS: Instrumental Parameters Development

C. A. Harvey, R. Meissner

November 8, 2005

## Disclaimer

---

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

This work was performed under the auspices of the U.S. Department of Energy by University of California, Lawrence Livermore National Laboratory under Contract W-7405-Eng-48.

## Pentaerythritol Tetranitrate (PETN) Surveillance by HPLC-MS: Instrumental Parameters Development

### Summary:

Surveillance of PETN Homologs in the stockpile here at LLNL is currently carried out by high performance liquid chromatography (HPLC) with ultra violet (UV) detection. Identification of unknown chromatographic peaks with this detection scheme is severely limited. The design agency is aware of the limitations of this methodology and ordered this study to develop instrumental parameters for the use of a currently owned mass spectrometer (MS) as the detection system. The resulting procedure would be a “drop-in” replacement for the current surveillance method (ERD04-524). The addition of quadrupole mass spectrometry provides qualitative identification of PETN and its homologs (Petrin, DiPEHN, TriPEON, and TetraPEDN) using a LLNL generated database, while providing mass clues to the identity of unknown chromatographic peaks.

### Experimental:

#### *Standards Preparation*

Stock standards were made in acetonitrile (ACN) from the following bulk materials.

- Pentaerythritol Tetranitrate (PETN): Lot B-331
- Pentaerythritol Trinitrate (PETRIN): Unknown lot, Russ Sanborn owner
- Dipentaerythritol Hexanitrate (DIPEHN): Unknown lot, Mound purified
- Tripentaerythritol Octanitrate (TRIPEON): Unknown lot, Mound purified
- Tetrapentaerythritol Decanitate (TETRAPEDN): Unknown lot, Mound purified

Calibration standards (50% ACN: 50% Water) were produced from these stocks in a concentration range from 0.5 mg/L to 22 mg/L

#### *Calibration*

Peak areas were recorded for both UV and MS detectors. Acceptable linearity was demonstrated with UV detection. The MS detector, while providing powerful qualitative information, did not provide linear response over the selected calibration range.

#### *HPLC Conditions*

After systematic evaluation of chromatographic parameters, the HPLC conditions were set as follows and are detailed in Appendix A.

- Column: Agilent LiChrospher RP-8, 125mm x 4mm ID, 100 Å pore size, 5 µm particle size.
- Eluent: 65% Acetonitrile and 35% Organic free Water with 0.05M Ammonium Formate as a modifier.
- Flow Rate: 0.25 ml/min.
- Column Temperature: 40°C.
- UV Detector Monitor: 204 nm, 4 nm bandwidth, 360 nm reference, 100nm bandwidth.
- Injection volume: 25 µl.

*Mass Spectrometer Conditions*

Run conditions of the LC/MSD are as follows. Additional detail is found in Appendix A.

- Ion polarity: Negative
- Ionization mode: Atmospheric Pressure Ionization-Electrospray (API-ES).
- Scan mass range: 50 to 1000 amu.
- Spray chamber
  - Gas Temperature: 250°C.
  - Drying Gas Flow: 10 L/min
  - Nebulizer Pressure: 50 psig
  - Capillary Voltage: 3000 V positive, 4000 V negative.

*Method Detection Limits (MDL)*

Standards representing sample homolog concentrations of 0.1 to 0.2 % were run seven times to determine the MDL for each analyte. The results were calculated by the following formula and are found in Table 1. Raw data in the form of chromatograms can be found in Appendix B.

$$MDL_{(7)} = s \cdot 3.14$$

s = Standard Deviation

Table 1: Method Detection Limits

<b>UV detector</b>	<b>PETRIN</b>	<b>PETN</b>	<b>DIPEHN</b>	<b>TRIPEON</b>	<b>TETRAPEDN</b>
Run 1	0.206159	0.120998	0.107454	0.107503	0.109236
Run 2	0.20391	0.112308	0.106364	0.104392	0.111241
Run 3	0.208033	0.118146	0.107324	0.10715	0.108766
Run 4	0.203682	0.119003	0.106294	0.105741	0.106282
Run 5	0.206018	0.118078	0.106296	0.107912	0.104756
Run 6	0.206406	0.115106	0.111418	0.11308	0.104134
Run 7	0.213287	0.116912	0.111839	0.107123	0.105963
Average	0.206785	0.117222	0.108141	0.107557	0.107197
Std. Dev.	0.003237	0.002821	0.002433	0.002717	0.002604
MDL (%)	<b>0.010</b>	<b>0.009</b>	<b>0.008</b>	<b>0.009</b>	<b>0.008</b>
<b>MS detector</b>	<b>PETRIN</b>	<b>PETN</b>	<b>DIPEHN</b>	<b>TRIPEON</b>	<b>TETRAPEDN</b>
Run 1	0.172987	0.165535	0.128442	0.126056	0.123501
Run 2	0.167646	0.160089	0.120958	0.117041	0.107316
Run 3	0.174055	0.1672	0.125997	0.1202	0.137822
Run 4	0.162921	0.130496	0.116445	0.128098	0.117446
Run 5	0.162582	0.156809	0.114358	0.123708	0.109512
Run 6	0.168198	0.181352	0.114748	0.119172	0.101854
Run 7	0.156901	0.190846	0.119167	0.121275	0.112809
Average	0.16647	0.164618	0.120016	0.122221	0.115751
Std. Dev.	0.006105	0.019241	0.005491	0.003928	0.011986
MDL (%)	<b>0.019</b>	<b>0.060</b>	<b>0.017</b>	<b>0.012</b>	<b>0.038</b>

*Data Analysis*

Data reduction and results calculation was done by Agilent Chemstation software. Calibration curves and integration parameters are found in the detailed method listing in Appendix C.

## Results and Discussion:

### *Method Sensitivity*

The MS detection platform proved unable to hit the RM Spec detection limit of 0.01%. However, the UV detector meets this requirement for all analytes. Since the UV detector is non-destructive, it is recommended that it be used in series prior to the MS. This will meet current RM Spec criteria while providing qualitative analyte identification.

### *Key Method Constraints*

The use of currently owned technology due to the limited scope and budget of this work led to the following constraints.

- Limiting the analysis to HPLC-MS eliminated the possibility of exploring other techniques potentially better suited to detect PETN breakdown products.
- The silica cartridge based column employed precludes the use of very acidic or basic conditions, limiting our chromatographic flexibility.
- Quadrupole mass spectrometry has limited mass accuracy ( $\pm 1$  amu), which prevents positive identification of molecular formula.
- HPLC-MS detection with electrospray introduction leads to adduction of the eluted compound to co-sprayed ions due to conditions found in the mass spectrometer source. The effective use of “user generated” libraries is highly subject to run conditions. Stringent control of method parameters and chemical quality is absolutely necessary to maintain functionality.

### *Recommendations*

The focus of this work was on producing a method with the currently owned quadrupole detection platform and chromatography system. It must be emphasized that this instrumentation is not the optimum set-up to generate the best liquid chromatography-mass spectrometry data available today. A broadened scope and funding is required to produce a truly “enhanced” methodology for the characterization of PETN powders. In my opinion, while somewhat underwhelming in overall performance, the procedure produced in this work is an adequate replacement for the current HPLC-UV method used in detonator surveillance (ERD04-524).

Currently published work<sup>(1)</sup>, using the exact same instrumentation employed in this work, showed no evidence of newly identified breakdown products. Additionally, work published here at LLNL suggests that potential breakdown products may not be physically amenable to HPLC analysis<sup>(2)</sup>. This suggests that other methodologies must be pursued to potentially unlock degradation pathways not yet documented in stockpile relevant samples. Industrial partnerships with companies willing to run demonstration samples on “state of the art” instrumentation could provide us with vital information to solidify the path an improved understanding of PETN aging in the enduring stockpile.

## References

1. Brackett C., “Contaminants and Decomposition Products in Naturally Aged PETN”, Dissertation, University of Pacific, 2005.
2. Chambers D. M., “Perspectives on Pentaerythritol Tetranitrate (PETN) Decomposition”, UCRL-ID-148956, July 1, 2002.

# Appendix A

## HPLC Run Method

Method Information

Homolog Content of Pentaerythritol Tetranitrate (PETN) by High  
Performance Liquid Chromatography/Mass Spectrometry (HPLC/MS)--65%  
Acetonitrile, 35% Water w/0.05M Ammonium Formate

Run Time Checklist

Pre-Run Cmd/Macro: off

Data Acquisition: on

Standard Data Analysis: on

Customized Data Analysis: off

Save GLP Data: off

Post-Run Cmd/Macro: off

Save Method with Data: off

## =====

## 1100 Binary Pump 1

## Control

Column Flow : 0.250 ml/min  
Stoptime : 50.00 min  
Posttime : Off

## Solvents

Solvent A 1 : 35.0 % (0.05M NH4COOH in Water)  
Solvent B 1 : 65.0 % (Acetonitrile)

## PressureLimits

Minimum Pressure : 10 bar  
Maximum Pressure : 300 bar

## Auxiliary

Maximal Flow Ramp : 100.00 ml/min^2  
Compressibility A : 50\*10^-6/bar  
Minimal Stroke A : Auto  
Compressibility B : 115\*10^-6/bar  
Minimal Stroke B : Auto

## Store Parameters

Store Ratio A : Yes  
Store Ratio B : Yes  
Store Flow : Yes  
Store Pressure : Yes

Timetable is empty

## =====

## Agilent 1100 Diode Array Detector 1

## Signals

Signal	Store	Signal, Bw	Reference, Bw	[nm]
A:	Yes	204 4	360 100	
B:	No	254 16	360 100	
C:	No	210 8	360 100	
D:	No	230 16	360 100	
E:	No	280 16	360 100	

## Spectrum

Store Spectra : Apex + Slope + Baselines  
Range from : 190 nm  
Range to : 400 nm  
Range step : 2.00 nm  
Threshold : 1.00 mAU

## Time

Stoptime : As pump  
Posttime : Off

## Required Lamps

UV lamp required : Yes  
Vis lamp required : No

## Autobalance

Prerun balancing : Yes  
Postrun balancing : No



Margin for negative Absorbance: 100 mAU

Peakwidth : > 0.1 min  
Slit : 4 nm

## Analog Outputs

Zero offset ana. out. 1: 5 %  
Zero offset ana. out. 2: 5 %  
Attenuation ana. out. 1: 1000 mAU  
Attenuation ana. out. 2: 1000 mAU

Timetable is empty

=====  
Mass Spectrometer Detector  
=====General Information  
-----

Use MSD : Enabled  
Ionization Mode : API-ES  
Tune File : atunes.tun  
StopTime : asPump  
Time Filter : Enabled  
Data Storage : Condensed  
Peakwidth : 0.10 min  
Fast Scan : Enabled  
Fast Scan Data Reconstruction: Enabled

Signals  
-----

[Signal 1]

Polarity : Negative  
Fragmentor Ramp : Disabled

## Scan Parameters

Time (min)	Mass Range Low High	Frag- mentor	Gain EMV	Thres- hold	Step- size
0.00	50.00 1000.00	40	4.0	150	0.10

[Signal 2]

Not Active

[Signal 3]

Not Active

[Signal 4]

Not Active

Spray Chamber  
-----

## [MSZones]

Gas Temp : 250 C  
DryingGas : 10.0 l/min  
Neb Pres : 50 psig  
Quad Temp : 0 C

VCap (Positive) : 3000 V  
VCap (Negative) : 4000 V

## [Time Table]

Time Table is empty.

UCRL-XX-XXXX  
Chemistry & Chemical Engineering Division  
November 4, 2005

## END OF MS ACQUISITION PARAMETERS

## =====

## FIA Series

FIA Series in this Method : Disabled

## Time Setting

Time between Injections : 2.00 min  
Injection Loop Flush Time : 0.17 min

## =====

## Agilent 1100 Autosampler 1

## Injection

Injection Mode : Standard  
Injector volume : 25.00 µl  
Optimization : none

## Auxiliary

Drawspeed : 100 µl/min  
Ejectspeed : 100 µl/min  
Draw position : 0.0 mm

## Time

Stoptime : As Pump  
Posttime : Off

## =====

## Agilent 1100 Column Thermostat 1

## Temperature settings

Left temperature : 40.0°C  
Right temperature : Not controlled  
Enable analysis : When Temp. is within setpoint +/- 0.8°C  
Store left temperature : Yes  
Store right temperature : Yes

## Time

Stoptime : As pump  
Posttime : Off

Column Switching Valve : Column 1

Timetable is empty



# Appendix B

## MDL Raw Data

Injection Date : 5/17/2005 6:44:29 PM

Seq. Line : 3

Sample Name : 1mg/l PETN+HL

Chemistry &amp; Chemical Engineering Division

Acq. Operator : R. Meissner

Inj : 1

November 4, 2005

Acq. Instrument : Instrument 1

Inj Volume : 25 µl

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M


Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)

Homolog Content of PETN by High Performance Liquid HPLC/MS



## External Standard Report

Sorted By : Signal  
 Calib. Data Modified : 5/6/2005 12:06:48 PM  
 Multiplier : 1.0000  
 Dilution : 1.0000  
 Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

RetTime [min]	Type	Area [mAU*s]	Amt/Area	Amount [ng/ul]	Grp	Name
7.404	VP	88.54317	2.32834e-2	2.06159		PETRIN
10.201	VB	207.41165	5.83369e-3	1.20998		PETN
15.746	VB	176.72702	6.08022e-3	1.07454		DIPEHN
25.024	VB	173.69185	6.18931e-3	1.07503		TRIPEON
39.807	VB	142.57104	7.66188e-3	1.09236		TETRAPEDN

Totals : 6.51350

Signal 2: MSD1 TIC, MS File

Injection Date : 5/17/2005 6:44:29 PM  
Sample Name : 1mg/l PETN+HL  
Acq. Operator : R. Meissner  
Acq. Instrument : Instrument 1  
Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M  
Last changed : 5/11/2005 1:35:14 PM by R. Meissner  
Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M  
Last changed : 5/18/2005 8:55:03 AM by R. Meissner  
(modified after loading)

Seq. Line : 3  
Location : Vial 3  
Inj : 1  
Inj Volume : 25 µl

UCRL-XX-XXXX  
Chemistry & Chemical Engineering Division  
November 4, 2005

Homolog Content of PETN by High Performance Liquid HPLC/MS

RetTime [min]	Type	Area	Amt/Area	Amount [ng/ul]	Grp	Name
7.556	VB	5.38170e6	3.52073e-7	1.89475		PETRIN
10.363	BP	2.94514e6	6.44215e-7	1.89730		PETN
15.867	VP	2.23655e7	5.67580e-8	1.26942		DIPEHN
25.240	VB	2.41875e7	5.29007e-8	1.27953		TRIPEON
39.947	BB	8.88557e6	1.39766e-7	1.24190		TETRAPEDN

Totals : 7.58291

\*\*\* End of Report \*\*\*

Injection Date : 5/17/2005 7:36:39 PM

Sample Name : 1mg/l PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)


Homolog Content of PETN by High Performance Liquid HPLC/MS

Seq. Line : 3  
Location : Vial 5Inj : 2  
Inj Volume : 25 µl

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005



## External Standard Report

Sorted By : Signal  
 Calib. Data Modified : 5/6/2005 12:06:48 PM  
 Multiplier : 1.0000  
 Dilution : 1.0000  
 Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

RetTime [min]	Type	Area [mAU*s]	Amt/Area	Amount [ng/ul]	Grp	Name
7.400	VP	87.57712	2.32834e-2	2.03910		PETRIN
10.193	VB	192.51671	5.83369e-3	1.12308		PETN
15.730	VB	174.93390	6.08022e-3	1.06364		DIPEHN
25.004	BP	168.66493	6.18931e-3	1.04392		TRIBEON
39.798	PB	145.18730	7.66188e-3	1.11241		TETRAPEDN

Totals : 6.38214

Signal 2: MSD1 TIC, MS File

Injection Date : 5/17/2005 7:36:39 PM  
Sample Name : 1mg/l PETN+HL  
Acq. Operator : R. Meissner  
Acq. Instrument : Instrument 1  
Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M  
Last changed : 5/11/2005 1:35:14 PM by R. Meissner  
Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M  
Last changed : 5/18/2005 8:55:03 AM by R. Meissner  
(modified after loading)

Seq. Line : 3  
Location : Vial 3  
Inj : 2  
Inj Volume : 25 µl

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005

Homolog Content of PETN by High Performance Liquid HPLC/MS

RetTime [min]	Type	Area	Amt/Area	Amount [ng/ul]	Grp	Name
7.537	PB	5.20483e6	3.52073e-7	1.83248		PETRIN
10.250	BP	3.28198e6	6.44215e-7	2.11430		PETN
15.839	VB	2.17218e7	5.67580e-8	1.23289		DIPEHN
25.067	BV	2.28299e7	5.29007e-8	1.20772		TRIPEON
39.931	BB	1.01905e7	1.39766e-7	1.42428		TETRAPEDN

Totals : 7.81167


\*\*\* End of Report \*\*\*



Injection Date : 5/17/2005 8:28:47 PM  
 Sample Name : 1mg/l PETN+HL  
 Acq. Operator : R. Meissner  
 Acq. Instrument : Instrument 1  
 Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M  
 Last changed : 5/11/2005 1:35:14 PM by R. Meissner  
 Analysis Method : C:\HPCHEM\1\METHODS\PETNROL.M  
 Last changed : 5/18/2005 8:55:03 AM by R. Meissner  
 (modified after loading)

Seq. Line : 3  
 Location : UCRL-XX-XXXX  
 Chemistry & Chemical Engineering Division  
 Inj : 3  
 Inj Volume : 25 µl  
 November 4, 2005

Homolog Content of PETN by High Performance Liquid HPLC/MS



### External Standard Report

Sorted By : Signal  
 Calib. Data Modified : 5/6/2005 12:06:48 PM  
 Multiplier : 1.0000  
 Dilution : 1.0000  
 Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

RetTime [min]	Type	Area [mAU*s]	Amt/Area	Amount [ng/ul]	Grp	Name
7.400	BV	89.34819	2.32834e-2	2.08033		PETRIN
10.198	VB	202.52371	5.83369e-3	1.18146		PETN
15.752	VB	176.51363	6.08022e-3	1.07324		DIPEHN
25.059	VB	173.12132	6.18931e-3	1.07150		TRIPEON
39.904	BP	141.95709	7.66188e-3	1.08766		TETRAPEDN

Totals : 6.49419

Signal 2: MSD1 TIC, MS File

Injection Date : 5/17/2005 8:28:47 PM

Sample Name : 1mg/1 PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)

Homolog Content of PETN by High Performance Liquid HPLC/MS

RetTime [min]	Type	Area	Amt/Area	Amount [ng/ul]	Grp	Name
7.544	BB	5.27285e6	3.52073e-7	1.85643		PETRIN
10.326	VP	3.38568e6	6.44215e-7	2.18111		PETN
15.844	VB	2.25486e7	5.67580e-8	1.27982		DIPEHN
25.111	BB	2.30834e7	5.29007e-8	1.22113		TRIPEON
39.887	BB	1.06411e7	1.39766e-7	1.48726		TETRAPEDN

Totals : 8.02573

1 Warnings or Errors :

Warning : Elution order of calibrated compounds may have changed

\*\*\* End of Report \*\*\*

Injection Date : 5/17/2005 9:21:01 PM

Sample Name : 1mg/l PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner


Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)

Homolog Content of PETN by High Performance Liquid HPLC/MS

Chemistry & Chemical Engineering Division  
 UCRL-XX-XXXX  
 November 4, 2005



## External Standard Report

Sorted By : Signal  
 Calib. Data Modified : 5/6/2005 12:06:48 PM  
 Multiplier : 1.0000  
 Dilution : 1.0000  
 Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

RetTime [min]	Type	Area [mAU*s]	Amt/Area	Amount [ng/ul]	Grp	Name
7.400	BP	87.47935	2.32834e-2	2.03682		PETRIN
10.196	VV	203.99232	5.83369e-3	1.19003		PETN
15.743	VB	174.81976	6.08022e-3	1.06294		DIPEHN
25.038	VB	170.84465	6.18931e-3	1.05741		TRIPLEON
39.870	BP	138.71533	7.66188e-3	1.06282		TETRAPEDN

Totals : 6.41002

Signal 2: MSD1 TIC, MS File

Injection Date : 5/17/2005 9:21:01 PM  
Sample Name : 1mg/1 PETN+HL  
Acq. Operator : R. Meissner  
Acq. Instrument : Instrument 1  
Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M  
Last changed : 5/11/2005 1:35:14 PM by R. Meissner  
Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M  
Last changed : 5/18/2005 8:55:03 AM by R. Meissner  
(modified after loading)

Seq. Line : 3  
Location : Vial 3  
Inj : 4  
Inj Volume : 25 µl

**UCRL-XX-XXXX**  
**Chemistry & Chemical Engineering Division**  
**November 4, 2005**

Homolog Content of PETN by High Performance Liquid HPLC/MS

RetTime [min]	Type	Area	Amt/Area	Amount [ng/ul]	Grp	Name
7.525	BB	5.64236e6	3.52073e-7	1.98652		PETRIN
10.274	BP	4.83963e6	6.44215e-7	3.11777		PETN
15.870	VB	2.28228e7	5.67580e-8	1.29538		DIPEHN
25.244	PV	2.51970e7	5.29007e-8	1.33294		TRIPEON
40.051	VB	9.37737e6	1.39766e-7	1.31064		TETRAPEDN

Totals : 9.04324

\*\*\* End of Report \*\*\*

Injection Date : 5/17/2005 10:13:08 PM

Sample Name : 1mg/l PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)

Homolog Content of PETN by High Performance Liquid HPLC/MS

Seq. Line : 3

Location : Vial 3


Inj : 5

Inj Volume : 25 µl

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005



## External Standard Report

Sorted By : Signal  
 Calib. Data Modified : 5/6/2005 12:06:48 PM  
 Multiplier : 1.0000  
 Dilution : 1.0000  
 Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

RetTime [min]	Type	Area [mAU*s]	Amt/Area	Amount [ng/ul]	Grp	Name
7.397	VP	88.48266	2.32834e-2	2.06018		PETRIN
10.194	VV	202.40721	5.83369e-3	1.18078		PETN
15.735	VP	174.82202	6.08022e-3	1.06296		DIPEHN
25.015	BB	174.35257	6.18931e-3	1.07912		TRIPLEON
39.802	BB	136.72397	7.66188e-3	1.04756		TETRAPEDN

Totals : 6.43060

Signal 2: MSD1 TIC, MS File

Injection Date : 5/17/2005 10:13:08 PM

Sample Name : 1mg/1 PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)

Homolog Content of PETN by High Performance Liquid HPLC/MS

Seq. Line : 3  
Location : Vial 3  
Inj : 5  
Inj Volume : 25 µl

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005

RetTime [min]	Type	Area	Amt/Area	Amount [ng/ul]	Grp	Name
7.517	PB	5.07498e6	3.52073e-7	1.78676		PETRIN
10.351	BP	3.63243e6	6.44215e-7	2.34007		PETN
15.856	VB	2.21522e7	5.67580e-8	1.25732		DIPEHN
25.157	BB	2.29474e7	5.29007e-8	1.21394		TRIPEON
40.628	VB	2.02553e6	1.39766e-7	2.83099e-1		TETRAPEDN

Totals : 6.88118

\*\*\* End of Report \*\*\*

Injection Date : 5/17/2005 11:05:20 PM

Sample Name : 1mg/1 PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)


Homolog Content of PETN by High Performance Liquid HPLC/MS

Seq. Line : 3  
Location : Vial 3  
Inj : 6  
Inj Volume : 25 µl

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005



## External Standard Report

Sorted By : Signal  
Calib. Data Modified : 5/6/2005 12:06:48 PM  
Multiplier : 1.0000  
Dilution : 1.0000  
Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

RetTime [min]	Type	Area [mAU*s]	Amt/Area	Amount [ng/ul]	Grp	Name
7.393	VP	88.64935	2.32834e-2	2.06406		PETRIN
10.184	VV	197.31152	5.83369e-3	1.15106		PETN
15.714	VP	183.24591	6.08022e-3	1.11418		DIPEHN
24.960	VV	182.70198	6.18931e-3	1.13080		TRIPLEON
39.578	VV	135.91180	7.66188e-3	1.04134		TETRAPEDN

Totals : 6.50143

Signal 2: MSD1 TIC, MS File

Injection Date : 5/17/2005 11:05:20 PM

Sample Name : 1mg/1 PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)

Homolog Content of PETN by High Performance Liquid HPLC/MS

Seq. Line : 3  
Location : Vial 3  
Inj : 6  
Inj Volume : 25 µl

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005

RetTime [min]	Type	Area	Amt/Area	Amount [ng/ul]	Grp	Name
7.542	PP	4.75882e6	3.52073e-7	1.67545		PETRIN
10.296	BP	3.06897e6	6.44215e-7	1.97708		PETN
15.798	BB	2.04323e7	5.67580e-8	1.15970		DIPEHN
25.176	BB	2.26709e7	5.29007e-8	1.19931		TRIPEON
39.600	BB	8.52627e6	1.39766e-7	1.19168		TETRAPEDN

Totals : 7.20321

\*\*\* End of Report \*\*\*



Injection Date : 5/17/2005 11:57:33 PM

Sample Name : 1mg/1 PETN+HL

Acq. Operator : R. Meissner

Acq. Instrument : Instrument 1

Acq. Method : C:\HPCHEM\1\METHODS\PETNLCMS.M

Last changed : 5/11/2005 1:35:14 PM by R. Meissner

Analysis Method : C:\HPCHEM\1\METHODS\PETNRO1.M

Last changed : 5/18/2005 8:55:03 AM by R. Meissner

(modified after loading)



Homolog Content of PETN by High Performance Liquid HPLC/MS

Seq. Line : 3  
 Location : Vial 3  
 Inj : 7  
 Inj Volume : 25 µl

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005



## External Standard Report

Sorted By : Signal  
 Calib. Data Modified : 5/6/2005 12:06:48 PM  
 Multiplier : 1.0000  
 Dilution : 1.0000  
 Use Multiplier & Dilution Factor with ISTDs

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

RetTime [min]	Type	Area [mAU*s]	Amt/Area	Amount [ng/ul]	Grp	Name
7.388	VP	91.60476	2.32834e-2	2.13287		PETRIN
10.172	VV	200.40895	5.83369e-3	1.16912		PETN
15.686	VB	183.93938	6.08022e-3	1.11839		DIPEHN
24.901	BB	173.07683	6.18931e-3	1.07123		TRIEPON
39.556	BB	138.29953	7.66188e-3	1.05963		TETRAPEDN

Totals : 6.55125

Signal 2: MSD1 TIC, MS File

=====

Injection Date	: 5/17/2005 11:57:33 PM	Seq. Line	: 3	UCRL-XX-XXXX
Sample Name	: 1mg/1 PETN+HL	Section	: 1	Chemistry & Chemical Engineering Division
Acq. Operator	: R. Meissner	Inj	: 7	November 4, 2005
Acq. Instrument	: Instrument 1	Inj Volume	: 25 µl	
Acq. Method	: C:\HPCHEM\1\METHODS\PETNLCMS.M			
Last changed	: 5/11/2005 1:35:14 PM by R. Meissner			
Analysis Method	: C:\HPCHEM\1\METHODS\PETNR01.M			
Last changed	: 5/18/2005 8:55:03 AM by R. Meissner			
	(modified after loading)			

Homolog Content of PETN by High Performance Liquid HPLC/MS

=====

RetTime [min]	Type	Area	Amt/Area	Amount [ng/ul]	Grp	Name
7.497	PB	4.40261e6	3.52073e-7	1.55004		PETRIN
10.287	BB	3.30767e6	6.44215e-7	2.13085		PETN
15.792	VB	2.08461e7	5.67580e-8	1.18318		DIPEHN
25.044	BB	2.34369e7	5.29007e-8	1.23983		TRIPEON
39.866	BB	7.84607e6	1.39766e-7	1.09661		TETRAPEDN

Totals : 7.20052

=====

\*\*\* End of Report \*\*\*

# Appendix C

## Data Reprocessing Method

Method Information

Homolog Content of PETN by High Performance Liquid HPLC/MS--Offline  
recalculation method.

## Estimated Sample Purity Calculation

November 4, 2005

## General

Estimated Sample Purity Calculation is disabled.

## Integration Events

## Default Integration Event Table "Event"

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

## Detector Default Integration Event Table "Event\_ADC"

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

## Detector Default Integration Event Table "Event\_FLD"

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

## Signal Specific Integration Event Table "Event\_MSD1SPC"

Event	Value	Time
Initial Slope Sensitivity(Full Scan)	1.000	Initial
Initial Peak Width(Full Scan)	0.250	Initial
Initial Slope Sensitivity(Cond. Scan/SIM)	0.100	Initial
Initial Peak Width(Cond. Scan/SIM)	0.050	Initial
Initial Area Reject	0.000	Initial
Initial Height Reject	5.000	Initial
Initial Shoulders	OFF	Initial

UCRL-XX-XXXX

Chemistry & Chemical Engineering Division  
November 4, 2005

Detector Default Integration Event Table "Event\_VWD"

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

Detector Default Integration Event Table "Event\_ECD"

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

Detector Default Integration Event Table "Event\_MWD"

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

Detector Default Integration Event Table "Event\_MSD"

Event	Value	Time
Initial Slope Sensitivity	10000.000	Initial
Initial Peak Width	0.070	Initial
Initial Area Reject	1000.000	Initial
Initial Height Reject	100.000	Initial
Initial Shoulders	OFF	Initial

Detector Default Integration Event Table "Event\_DAD"

Event	Value	Time
Initial Slope Sensitivity	0.515	Initial
Initial Peak Width	0.089	Initial
Initial Area Reject	10.000	Initial
Initial Height Reject	0.028	Initial
Initial Shoulders	OFF	Initial
Integration	OFF	0.000
Integration	ON	2.000

UCRL-XX-XXXX

Chemistry & Chemical Engineering Division  
November 4, 2005

## Signal Specific Integration Event Table "Event\_DAD11"

Event	Value	Time
Initial Slope Sensitivity	0.100	Initial
Initial Peak Width	0.200	Initial
Initial Area Reject	10.000	Initial
Initial Height Reject	0.028	Initial
Initial Shoulders	OFF	Initial
Integration	OFF	0.000
Integration	ON	6.000

## Signal Specific Integration Event Table "Event\_MSD1TIC"

Event	Value	Time
Initial Slope Sensitivity	90000.000	Initial
Initial Peak Width	0.500	Initial
Initial Area Reject	20000.000	Initial
Initial Height Reject	10000.000	Initial
Initial Shoulders	OFF	Initial
Integration	OFF	0.000
Integration	ON	6.000

Apply Manual Integration Events: No

Advanced Baseline : No

## Specify Report

Destination: Screen  
 Quantitative Results sorted by: Signal  
 Report Style: Short  
 Sample info on each page: Yes  
 Add Chromatogram Output: Yes  
 Chromatogram Output: Portrait  
 Size in Time direction: 100 % of Page  
 Size in Response direction: 40 % of Page

## Signal Options

Include: Axes, Retention Times, Baselines, Tick Marks  
 Font: Arial, Size: 8

Ranges: Use Ranges

	Min Value	Max Value
Time	0.000	50.000
Response		

Multi Chromatograms: Separated, Each in full Scale

## Calibration Table

Calib. Data Modified : 5/6/2005 12:06:48 PM

29 of 31

UCRL-XX-XXXX

Chemistry &amp; Chemical Engineering Division

November 4, 2005

Calculate : External Standard  
 Based on : Peak Area  
 Rel. Reference Window : 5.000 %  
 Abs. Reference Window : 0.000 min  
 Rel. Non-ref. Window : 5.000 %  
 Abs. Non-ref. Window : 0.000 min  
 Use Multiplier & Dilution Factor with ISTDs  
 Uncalibrated Peaks : not reported  
 Partial Calibration : Yes, identified peaks are recalibrated  
 Correct All Ret. Times: No, only for identified peaks

Curve Type : Average Response/Amount  
 Origin : Ignored  
 Weight : Equal

Recalibration Settings:  
 Average Response : Average all calibrations  
 Average Retention Time: Floating Average New 75%

## Calibration Report Options :

Printout of recalibrations within a sequence:

Calibration Table after Recalibration

Normal Report after Recalibration

If the sequence is done with bracketing:

Results of first cycle (ending previous bracket)

Signal 1: DAD1 A, Sig=204,4 Ref=360,100

Signal 2: MSD1 TIC, MS File

RetTime	Lvl	Amount	Area	Amt/Area	Ref Grp Name
[min]	Sig	[ng/ul]			
7.411	1 4	3.20200e-1	21.10326	1.51730e-2	PETRIN
	3	8.00500e-1	32.30188	2.47818e-2	
	1	1.60100	57.62590	2.77826e-2	
	2	3.20200	117.37147	2.72809e-2	
	5	8.00500	295.15897	2.71210e-2	
	10	16.01000	717.93854	2.23000e-2	
	20	32.02000	1281.44275	2.49875e-2	
7.569	2 4	3.20200e-1	1.71619e6	1.86576e-7	PETRIN
	3	8.00500e-1	2.67383e6	2.99383e-7	
	1	1.60100	3.95461e6	4.04844e-7	
	2	3.20200	8.55226e6	3.74404e-7	
	5	8.00500	1.75550e7	4.55996e-7	
	10	16.01000	3.41416e7	4.68929e-7	
	20	32.02000	5.49389e7	5.82829e-7	
10.244	1 4	2.12800e-1	39.26180	5.42003e-3	PETN
	3	5.32000e-1	92.09941	5.77637e-3	
	1	1.06400	164.72203	6.45937e-3	
	2	2.12800	336.47046	6.32448e-3	
	5	5.32000	852.20819	6.24261e-3	
	10	10.64000	2082.84595	5.10840e-3	
	20	21.28000	3690.59302	5.76601e-3	
10.292	2 4	2.12800e-1	2.56349e4	8.30118e-6	PETN
	3	5.32000e-1	6.47165e5	8.22047e-7	
	1	1.06400	2.00867e6	5.29704e-7	
	2	2.12800	4.47120e6	4.75935e-7	
	5	5.32000	1.13131e7	4.70250e-7	
	10	10.64000	2.01378e7	5.28359e-7	
	20	21.28000	3.23634e7	6.57533e-7	
15.867	1 4	2.01800e-1	38.17426	5.28629e-3	DIPEHN
	3	5.04500e-1	85.46818	5.90278e-3	
	1	1.00900	152.95752	6.59660e-3	
	2	2.01800	303.11578	6.65752e-3	
	5	5.04500	754.09277	6.69016e-3	
	10	10.09000	1835.50110	5.49714e-3	

30 of 31



						UCRL-XX-XXXX	
RetTime	Lvl	Amount	Area	Amt/Area	Ref Grp Name		
[min]	Sig	[ng/ul]			Chemistry & Chemical Engineering Division		
-----	--	----	-----	-----	----	--	-----
							November 4, 2005
15.969	2	20 20.18000	3218.86792	6.26928e-3			
		4 2.01800e-1	5.40204e6	3.73562e-8	DIPEHN		
		3 5.04500e-1	1.19974e7	4.20506e-8			
		1 1.00900	1.88990e7	5.33891e-8			
		2 2.01800	3.35915e7	6.00748e-8			
		5 5.04500	6.84058e7	7.37511e-8			
		10 10.09000	1.34841e8	7.48290e-8			
25.237	1	20 20.18000	2.11510e8	9.54092e-8			
		4 2.02000e-1	38.45569	5.25280e-3	TRIPEON		
		3 5.05000e-1	83.18657	6.07069e-3			
		1 1.01000	145.42653	6.94509e-3			
		2 2.02000	298.88165	6.75853e-3			
		5 5.05000	736.34021	6.85824e-3			
		10 10.10000	1813.55432	5.56917e-3			
		20 20.20000	3203.00122	6.30659e-3			
25.492	2	4 2.02000e-1	5.94654e6	3.39694e-8	TRIPEON		
		3 5.05000e-1	1.28702e7	3.92379e-8			
		1 1.01000	1.83797e7	5.49519e-8			
		2 2.02000	3.69940e7	5.46034e-8			
		5 5.05000	7.58113e7	6.66128e-8			
		10 10.10000	1.47407e8	6.85177e-8			
		20 20.20000	2.27876e8	8.86447e-8			
39.974	1	4 2.24200e-1	34.64077	6.47214e-3	TETRAPEDN		
		3 5.60500e-1	80.59784	6.95428e-3			
		1 1.12100	131.03842	8.55474e-3			
		2 2.24200	246.60544	9.09145e-3			
		5 5.60500	599.11951	9.35540e-3			
		10 11.21000	1639.77954	6.83628e-3			
		20 22.42000	3032.37524	7.39354e-3			
40.275	2	4 2.24200e-1	4.27958e5	5.23883e-7	TETRAPEDN		
		3 5.60500e-1	6.52629e6	8.58835e-8			
		1 1.12100	8.79406e6	1.27472e-7			
		2 2.24200	1.82329e7	1.22964e-7			
		5 5.60500	3.67238e7	1.52626e-7			
		10 11.21000	8.36064e7	1.34081e-7			
		20 22.42000	1.46712e8	1.52816e-7			

=====

Peak Sum Table

=====

\*\*\*No Entries in table\*\*\*

=====