



LAWRENCE
LIVERMORE
NATIONAL
LABORATORY

Diagnostic evaluation of a multiplexed RT-PCR microsphere array assay for the detection of foot-and-mouth disease virus and look-alike disease viruses

B. J. Hindson, S. M. Reid, B. R. Baker, K. Ebert, N. P. Ferris, L. F. Bentley Tammoro, R. J. Lenhoff, P. Naraghi-Arani, E. A. Vitalis, T. R. Slezak, P. J. Hullinger, D. P. King

August 1, 2007

Journal of Clinical Microbiology

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

**Diagnostic evaluation of a multiplexed RT-PCR microsphere array assay for the detection of foot-
and-mouth and look-alike disease viruses**

Running title: Multiplexed RT-PCR detection of FMDV and look-alikes

Benjamin J. Hindson^{1*}, Scott M. Reid², Brian R. Baker¹, Katja Ebert², Nigel P. Ferris², Lance F.
Bentley Tammemo¹, Raymond J. Lenhoff¹, Pejman Naraghi-Arani¹, Elizabeth A. Vitalis¹, Thomas R.
Slezak¹, Pamela J. Hullinger¹ and Donald P. King²

¹Lawrence Livermore National Laboratory, Livermore, California

²Institute for Animal Health, Pirbright Laboratory, Pirbright, Woking, Surrey, GU24 0NF, United
Kingdom

*Corresponding author. Mailing address: Lawrence Livermore National Laboratory, P.O. Box 808,
L-211, Livermore, CA 94551. Phone: (925) 423-8667. Fax: (925) 424-2778. E-mail: hindson1@llnl.gov

14 **ABSTRACT**

15 A high-throughput multiplexed assay was developed for the differential laboratory diagnosis of foot-
16 and-mouth disease virus (FMDV) from viruses which cause clinically similar diseases of livestock.
17 This assay simultaneously screens for five RNA and two DNA viruses using multiplexed reverse
18 transcription PCR (mRT-PCR) amplification coupled with a microsphere hybridization array and flow-
19 cytometric detection. Two of the seventeen primer-probe sets included in this multiplex assay were
20 adopted from previously characterized real-time RT-PCR (rRT-PCR) assays for FMDV. The diagnostic
21 accuracy of the mRT-PCR was evaluated using 287 field samples, including 248 (true positive n= 213,
22 true negative n=34) from suspect cases of foot-and-mouth disease collected from 65 countries between
23 1965 and 2006 and 39 true negative samples collected from healthy animals. The mRT-PCR assay
24 results were compared with two singleplex rRT-PCR assays, using virus isolation with antigen-ELISA
25 as the reference method. The diagnostic sensitivity of the mRT-PCR assay for FMDV was 93.9% [95%
26 C.I. 89.8-96.4%], compared to 98.1% [95% C.I. 95.3-99.3%] for the two singleplex rRT-PCR assays
27 used in combination. In addition, the assay could reliably differentiate between FMDV and other
28 vesicular viruses such as swine vesicular disease virus and vesicular exanthema of swine virus.
29 Interestingly, the mRT-PCR detected parapoxvirus (n=2) and bovine viral diarrhea virus (n=2) in
30 clinical samples, demonstrating the screening potential of this mRT-PCR assay to identify viruses in
31 FMDV-negative material not previously recognized using focused single-target rRT-PCR assays.

32 INTRODUCTION

33 Foot-and-mouth disease (FMD) is a highly infectious and contagious vesicular disease affecting
34 domestic and wild ruminants and swine caused by a single-stranded positive-sense RNA virus having
35 seven distinct serotypes (A, Asia 1, C, O, SAT 1, SAT 2 and SAT 3) (17). FMD is endemic in many
36 countries throughout the world, with serotype O having the highest prevalence, followed by serotype A
37 (9). Early detection of the virus is critical to minimizing disease spread and the significant economic
38 implications (12) resulting from the introduction of FMD into a country previously free of the disease.
39 Diagnosis of FMD can be confounded by diseases with similar clinical signs (“look-alike” diseases),
40 and by species where presentation of the disease is mild or indistinct (16). For the laboratory
41 identification of FMDV, the *Office International des Epizooties* (OIE) recommends virus isolation (VI),
42 antigen-ELISA (Ag-ELISA), and RT-PCR with detection by agarose gel electrophoresis or in real-time
43 using TaqMan® fluorogenic probes (15).

44 Real-time PCR is widely used by diagnostic laboratories, to complement or as a replacement for more
45 traditional detection methods. Two independent real-time reverse transcription polymerase chain
46 reaction (rRT-PCR) assays for FMD laboratory diagnosis target the ribosomal entry site of the
47 5′ untranslated region (5′UTR) (30) and the viral RNA polymerase gene (3D) (4) on the highly variable
48 FMDV genome. The 5′UTR and 3D rRT-PCR assays were initially compared prior to their
49 implementation in Australia (3). A subsequent in-depth comparative evaluation (14) was conducted to
50 further evaluate the effectiveness of these assays; demonstrating a higher diagnostic sensitivity of the
51 rRT-PCR assays over VI and/or antigen-ELISA (35), particularly when both assays were used in
52 combination. Both assays are used routinely in combination at the Food and Agriculture Organization
53 of the United Nations, World Reference Laboratory (FAO WRL) for FMD. rRT-PCR assays have also
54 been reported for the detection of other viruses which cause vesicular disease of livestock including
55 swine vesicular disease (SVD) (29), vesicular stomatitis (VS) (11, 28) and vesicular exanthema of swine
56 (VES) (31) or symptomatic look-alike diseases including bluetongue (13, 27, 33), bovine viral diarrhea

57 (1, 2, 20, 22, 39), malignant catarrhal fever (37) and parapoxvirus (26).

58 Fluorescent probes for rRT-PCR detection have broad emission spectra which limit multiplexing
59 capacity to the four or five discrete optical channels typically present in most commercial real-time PCR
60 instruments. Therefore, simultaneous testing for FMDV and look-alike disease viruses by rRT-PCR
61 would require many assays to be run in parallel, thereby increasing the demand on instrumentation and
62 reagents which escalates costs. A single multiplexed screening test that simultaneously detects and
63 differentiates FMDV from look-alike disease viruses would be desirable. Such a test could facilitate
64 rapid and cost-effective screening of suspect FMD field samples for laboratory differential diagnosis,
65 targeted FMD surveillance, or embedded foreign animal disease surveillance whilst conducting routine
66 testing for endemic diseases.

67 Luminex xMAP technology is a multiplexed high-throughput detection system (38) with many
68 applications for nucleic acid detection (6). The Luminex array offers up to 100 independent channels
69 using microspheres (5.6 μm) embedded with varying ratios of two fluorescent dyes. User-defined
70 surface modifications can include oligonucleotides, antibodies, peptides or other macromolecules.
71 Typically, a mixed suspension of functionalized microspheres is mixed with the sample to bind analytes
72 which are then labeled with a fluorescent reporter and analyzed using a specialized flow-cytometer. The
73 assay provides end-point detection with qualitative results by comparing fluorescence responses
74 (median fluorescence intensity, MFI) of each microsphere class to cut-off values. Recent nucleic acid
75 applications of the Luminex array include the detection and differentiation of Classical swine fever virus
76 from other pestiviruses (5), human respiratory viruses (18, 19, 21, 23), human papillomavirus (10, 32)
77 and human influenza A virus typing (40).

78 This report describes a novel multiplexed RT-PCR microsphere array assay for the differential
79 detection of FMDV from look-alike disease viruses. The development, optimization and analytical
80 evaluation of this multiplex assay will be reported separately. The primary purpose of this study was to
81 evaluate the diagnostic performance of the mRT-PCR assay for the detection of FMDV using a panel of

suspect field samples. The diagnostic sensitivities of the 3D and 5'UTR FMDV assays in mRT-PCR format were compared to rRT-PCR using VI with Ag-ELISA as the reference method. The detection of look-alike diseases in suspect FMD field samples is also reported.

MATERIALS AND METHODS

Samples. The panel comprised epithelia (true positive n= 213, true negative n=34) from suspect field cases of FMD submitted from 65 countries to the FAO WRL for FMD between 1965 and 2006 and included representatives of all seven serotypes of FMDV. In addition, 39 samples of true negative tongue epithelial were collected from healthy cattle at a UK abattoir. All sample testing and reference measurements were conducted at the FAO WRL for FMD. Epithelia were ground and suspended to generate a ~10% w/v suspension in phosphate buffer (0.04M, pH 7.6). The epithelial suspensions (ES) were centrifuged, the supernatant collected, then stored at -80°C. FMDV serotype was determined by conducting an antigen-ELISA (8) on the original ES or after viral propagation in cell culture. Samples of look-alike viruses from the FAO WRL for FMD collection included SVDV, VESV, San Miguel sea lion virus (SMSV), caliciviruses isolated from a variety of species (including cetacean, bovine, feline, reptile, skunk) and vesicular stomatitis virus (VSV).

Nucleic acid extraction. Total nucleic acid was extracted from each ES by an automated procedure using a MagNA Pure LC (Roche, UK) as previously described (14, 35). Extracted samples (40 µL) were aliquoted (3×13 µL), stored at -80°C and thawed once just before use.

rRT-PCR. Previously reported protocols for the individual 3D (4) and 5'UTR (34) rRT-PCR assays were modified for use in this study. Briefly, 25 µl reaction mixes (SuperScriptTM III Platinum[®] One-Step qRT-PCR System (Invitrogen)) containing 20 pmol of each primer, 7.5 pmol of dual-labelled TaqMan[®] probe and 5 µl total nucleic acid were prepared in an optical reaction plate (Stratagene, Amsterdam, The Netherlands). For both targets, RT-PCR amplification was performed in an Mx4000 Multiplex Quantitative PCR System (Stratagene) as described previously (34).

107 **mRT-PCR assay design.** A schematic depiction of the mRT-PCR assay is shown in Figure 1. The
108 RT-PCR uses eighteen biotinylated forward and unmodified reverse primer sets (17 for detection, 1
109 serves as a control). The sequences that comprise the multiplex assay are shown Table 1. The multiplex
110 assay was designed to detect and differentiate FMDV from SVDV, VESV, Bovine Viral Diarrhea Virus
111 (BVDV), Bluetongue Virus (BTV), Parapoxviruses (Orf virus, pseudocowpox virus, and bovine papular
112 stomatitis virus), and Bovine Herpes Virus-1 (BHV-1). Primer-probe sequences for FMDV (4, 30) and
113 BVDV (7, 22) were based on the work of others and adapted to the current multiplex format. All other
114 sequences were designed at Lawrence Livermore National Laboratory (LLNL) using an approach that
115 has previously been described (36). Additional computational analyses were performed to ensure
116 specificity and reliability against all available data, including a BLAST-based comparison of each
117 primer-probe set as a triplet against all sequences in GenBank to identify the targets that are predicted to
118 produce a PCR or TaqMan reaction at 57 °C for primer annealing and 67 °C for probe annealing, where
119 temperatures are derived from Primer 3 oligo T_M calculations. Optimal candidate primer-probe sets
120 were forwarded to the bench screening phase for further down-selection. Amplicon sizes ranged from
121 95 -349 bp.

122 **mRT-PCR primers and probes.** All oligonucleotides used for mRT-PCR were synthesized by
123 Integrated DNA Technologies (IDT, Coralville, IA) and HPLC-purified. Forward primers were
124 functionalized with 5' and internal biotin moieties. Reverse primers were unmodified. Probes were 5'
125 amino C6-modified with an 18-atom hexaethyleneglycol spacer. Lyophilized probe was dissolved in 2-
126 (*N*-morpholino)ethanesulfonic acid (MES) to yield a stock concentration of 1 mM. Lyophilized forward
127 and reverse primers were dissolved in TE buffer to yield a stock concentration of 1 mM. Working
128 dilutions were prepared from the stock solutions as required.

129 **Coupling of probe oligonucleotides to microspheres.** xMAP® Multi-Analyte COOH Microspheres
130 (Luminex Corp., Austin, TX) were covalently coupled to probe oligonucleotides using carbodiimide
131 activation based on the manufacturer's protocol. Briefly, stock microspheres (1 mL, 1.25×10^7

132 microspheres) were vortexed for 30 s, sonicated for 60 s, centrifuged at 8000 ×g for 5 min and the
133 supernatant removed. The microspheres were then resuspended in MES (50 µL, 0.1M, pH 4.5),
134 vortexed and sonicated. Probe (10 µL, 50 mM in MES) was added and the mixture vortexed. An
135 aqueous solution of *N*-(3-dimethylaminopropyl)-*N'*-ethylcarbodiimide hydrochloride (EDC; 5 µL, 10
136 mg/mL) was added, vortexed, then gently agitated for 30 min in the dark. A second aliquot of EDC (5
137 µL, 10 mg/mL) was added, vortexed, then gently agitated for 30 min in the dark. Tween 20 (1 mL,
138 0.02% v/v) was added, vortexed, centrifuged, the supernatant removed, then repeated using SDS (1 mL,
139 0.1% m/v) then TE buffer (100 µL 1 mM Tris, 10 mM EDTA, pH 7.4). The probe-conjugated
140 microspheres were resuspended in TE buffer (250 µL), vortexed, then stored at 4°C in the dark.

141 **Microsphere mixture.** A 21-plex microsphere suspension was prepared by combining individual
142 stock suspensions (9 µL/class) with Tris-NaCl buffer (3 mL). The mixture was vortexed then
143 enumerated using the Bio-Plex, targeting approximately 150 microspheres counts per class in 40 s. If
144 required, additional microspheres from the individual stocks were added to ensure the concentrations of
145 all classes were approximately equal.

146 **mRT-PCR amplification.** Each field sample was analyzed in duplicate by the multiplex assay.
147 Amplification was performed using a one-step RT-PCR kit (SuperScript™ III One-Step RT-PCR
148 System with Platinum® Taq DNA Polymerase, Invitrogen). The reaction volume of 25 µL was
149 comprised of nuclease-free water (0.95 µL), primer mix (3.6 µL), SuperScript™ III 2X reaction mix
150 (12.5 µL), MgSO₄ (0.95 µL, 50 mM, Invitrogen), SuperScript™ III RT / Platinum® Taq Mix (1 µL),
151 internal control armored RNA (1µL, ~100 copies) and template (5 µL). The internal control armored
152 RNA was lysed by heating to 70°C for 4 min prior to addition to the mix. The final concentration of
153 each primer and MgSO₄ was 0.4 µM and 3.5 mM, respectively. The RT-PCR thermal cycling protocol
154 was 55°C for 30 min, 95°C for 2 min, then 35 cycles of 95°C for 15 s, 60°C for 30 s, and 72°C for 15s,
155 followed by 72°C for 2 min with a final 4°C hold.

156 **Microsphere array hybridization.** A wash-assay format was adopted to reduce variability of the
157 response caused by non-specific hybridization of PCR products and fluorescent label. In a 96-well
158 plate, RT-PCR product (1 μ L) was added to a mixed suspension of probe-conjugated microspheres (22
159 μ L) then placed in a thermal cycler and subjected to 95°C for 2 min, 55°C for 5 min, followed by a 4°C
160 hold. Tris-NaCl buffer (100 μ L, 0.1M Tris, 0.2 M NaCl, 0.05% v/v Triton X-100, pH 8.0, Teknova) was
161 added and the suspension was transferred to a 96-well vacuum filter plate (MABVN 1250 Multiscreen
162 Filter Plate, Millipore). The suspension was vacuum-aspirated, and washed twice with Tris-NaCl buffer
163 (2 \times 100 μ L). Stock streptavidin phycoerythrin (SAPE, 1 mg/mL, Caltag Laboratories) was diluted with
164 Tris-NaCl to a working concentration of 3 μ g/mL. SAPE (60 μ L, 3 μ g/mL) was added then incubated
165 in the dark for 5 min. The suspension was vacuum-aspirated, washed once with Tris-NaCl (100 μ L), re-
166 suspended in Tris-NaCl (100 μ L) then transferred to a 96-well round bottom plate for fluorescent
167 detection. Nucleic acid extraction, PCR reaction assembly and PCR amplification were conducted in
168 separate rooms. Likewise, hybridization and Bio-Plex detection were performed together in a separate
169 room to minimize the likelihood of PCR contamination by amplicons.

170 **mRT-PCR detection.** Fluorescence detection of the processed microsphere suspension array was
171 achieved using a Bio-Plex Workstation (Bio-Rad, CA) set to count a minimum of 100 events per
172 microsphere class in a 50 μ L Bio-Plex sample volume. The Bio-Plex Workstation is a specialized dual-
173 laser flow cytometer integrated with XY microplate platform configured to analyze Luminex xMAP®
174 microspheres in a 96-well plate format. The reporter PMT voltage was calibrated on the low setting
175 with background subtraction enabled. The Bio-Plex workstation was validated and calibrated per the
176 manufacturer's instructions. The resolution of the MFI was 0.5 units for all channels.

177 **mRT-PCR assay controls.** The multiplex assay incorporates four control channels integral to each
178 reaction that are used to verify assay integrity. The negative control (NC) is a microsphere conjugated
179 to *Thermotoga maritima*-derived oligonucleotide sequence (MT-7) that serves as a non-specific binding
180 control in the multiplex PCR assay, and its response should remain consistently low (MFI \leq 80). The

181 fluorescence control (FC), biotinylated MT-7, confirms that fluorescent labeling with SAPE occurred
182 and should exhibit a high response (MFI >1000). The instrument control (IC) comprises a Cy3-labeled
183 MT-7 conjugate. As Cy3 and SAPE have similar fluorescence excitation and emission wavelengths, the
184 IC confirms proper function of the reporter optics within the Bio-Plex flow cytometer (MFI>500).
185 Armored RNA served as an end-to-end amplification control (AC) to reduce the probability of false
186 negative and is utilized at low concentrations (100 copies/reaction) to generate a low-level response
187 (MFI \geq 20) which minimizes its competition with detection channels. The armored RNA (XenoRNA-
188 01, Ambion, Austin, TX) is a proprietary 1070 nucleotide RNA transcript consisting of unique
189 nucleotide sequences that possess no significant homology to the current annotated sequences in
190 commonly used sequence databases including NCBI, Affymetrix, and Rosetta. Primers and probe-
191 labeled microsphere for the AC are included in the multiplex primer mix and microsphere suspension,
192 respectively.

193 **mRT-PCR assay integrity.** Minimum bead count quota and control channel responses were used to
194 verify the integrity of the detection channel responses which minimizes the likelihood of false positives
195 and false-negative results caused by operator error, instrument malfunction, non-specific hybridization,
196 or PCR inhibition. The responses of all detection channels are checked against pre-established cut-offs
197 then identified as mRT-PCR positive, negative or inconclusive. For each sample, the MFI of each
198 control bead class was checked against a cut-off value. A given result was considered invalid when the
199 MFI value of the IC, NC or FC controls was out of range, or when both the response of the AC was <20
200 and no detection channels exceeded cut-off. Results were considered valid when the AC was <20 and
201 any detection channel exceeded cut-off. The AC response can be diminished by a strong positive
202 response on a detection channel caused by competition in the RT-PCR reaction. If the responses of all
203 control channels are acceptable, the number of beads counted for control and detection channels are
204 checked (\geq 40 beads per channel). If a low bead count for any control channel occurred, then all results
205 for that sample were considered invalid. A low bead count on a given detection channel was considered

206 an invalid result for the channel in question.

207 **Data analysis.** Raw data exported from the Bio-Plex instrument were imported into MATLAB
208 (MathWorks) then analyzed using Microsoft Excel. If at least one result from a duplicate sample
209 analysis exceed cut-off, the sample was assigned as mRT-PCR positive. Receiver Operating
210 Characteristic (ROC) plots (24) were generated using a custom MATLAB program. Published
211 guidelines (25) were followed for the calculations of diagnostic test accuracy and statistical methods to
212 quantify uncertainty.

213 **Cut-off values.** The mRT-PCR responses for this assay are typically non-Gaussian and therefore
214 non-parametric methods were used to determine cut-offs. Each channel in the multiplex has a distinct
215 distribution of responses to true negative samples and therefore each is assigned its own cut-off value.
216 For each detection channel of the mRT-PCR assay, the responses to true negative samples were ranked
217 according to magnitude, then cut-off values were identified as the response (MFI value) which gave a
218 false positive rate closest to 5% (without exceeding 5%), corresponding to a diagnostic specificity of at
219 least 95% (3D MFI ≥ 6.5 , 5'UTR MFI ≥ 5.5). Other mRT-PCR detection channels for the FMD look-
220 alike disease viruses had similar cut-off values. For rRT-PCR, the cut-off was Ct ≤ 32 for both 3D and
221 5'UTR assays (34).

222

223 **RESULTS**

224

225 **Clinical sample validation.** Receiver Operating Characteristic (ROC) plots were constructed to
226 compare the diagnostic performance of the 3D and 5'UTR assays in rRT-PCR and mRT-PCR formats
227 (Figure 2). For the purpose of this evaluation, VI with Ag-ELISA served as the reference method and
228 therefore assigned samples as true negative or true positive, with the caveat that rRT-PCR is known to
229 detect FMDV in some samples considered negative by VI with Antigen-ELISA (35). Each plot shows
230 the true positive fraction (TPF; sensitivity) versus the false positive fraction (FPF; 1-specificity) over the

entire range of cut-off values. The ROC plots indicate that the 3D and 5'UTR assays in mRT-PCR format lost some ability to distinguish between true negative and true positive samples, as compared to rRT-PCR. The associated areas under each ROC curve also indicate the level of test performance in the absence of a cut-off value. Transfer to the multiplex format had the greatest effect upon the 5'UTR assay where the area under the curve was reduced to 0.773 in comparison to 0.942 for the rRT-PCR format. In contrast, the effect upon the 3D assay was less apparent: 0.955 and 0.985 for the mRT-PCR and rRT-PCR formats respectively. The ROC plots were generated using 3D and 5'UTR assay results independently. In practice, 3D and 5'UTR assay results would be used in combination, an intrinsic feature of the mRT-PCR format. The trade-off between TPF and FPF shown by ROC plots can be used to inform cut-off selection. The cut-off value for all detection channels was determined from the mRT-PCR response to true negative samples (n=74) using a specificity of 95%.

Table 2 summarizes the performance metrics of the 3D and 5'UTR assays when used independently or in combination for both formats. Results are presented according to serotype then summarized for all serotypes. In some cases, the rRT-PCR response for true positive samples (3D; n=5, 5'UTR; n=9) and true negative samples (3D; n=5, 5'UTR; n=2) yielded Ct values that were beyond the cut-off. In practice, these samples would be considered weak rRT-PCR positives and retested.

Table 3 shows a three-way comparison of mRT-PCR and rRT-PCR against the reference method using the combined 3D and 5'UTR assay results. The agreement between mRT-PCR and rRT-PCR for true positive samples was 95.8% (204/213). Two samples classified as FMDV-negative by VI and Ag-ELISA tested positive by both formats (SYR 6/2002, LAO 16/2003). Similarly, three samples (SYR 7/2002, TUR 17/2002, BHU 5/2004) which were negative by VI and Ag-ELISA tested positive by rRT-PCR but were negative by mRT-PCR due to its higher LOD. These findings are consistent with earlier studies (14, 35) where higher analytical sensitivity of rRT-PCR enabled detection of FMDV in samples designated negative by VI and Ag-ELISA. Due to its higher LOD, the mRT-PCR missed 9/209 positive samples detected by rRT-PCR. Four mRT-PCR false positives had MFI responses that were close to the

cut-off which was defined using a specificity of 95%. Four mRT-PCR and rRT-PCR false negatives were samples of FMDV serotype A (NIG 12/74), C (PHI 2/89) and O (O₁ Manisa TUR 8/69, YEM 1/2001). However, the rRT-PCR did generate Ct responses for two of these samples including PHI 2/89 (3D; 38.94, 5'UTR; 38.51) and O₁ Manisa (3D; 32.71). As the Ct values were beyond the cut-off they would be considered weak rRT-PCR positives and retested. These false negatives which generated either weak or no responses in mRT-PCR and rRT-PCR formats were likely caused by mutations in the FMDV gene segments targeted by the 3D and 5'UTR assays. They could also be due to low amounts of virus present, particularly for O₁ Manisa, which is a titrated sample used as a positive control in the diagnostic rRT-PCR that was deliberately set to be only weakly positive.

In addition to the 287 samples used to evaluate the mRT-PCR assay, eleven additional suspect FMDV field samples of undetermined serotype were analyzed. These samples were found to be negative by VI and Ag-ELISA but confirmed positive by combined 3D and 5'UTR rRT-PCR during a previous study. These samples were analyzed by the mRT-PCR assay, which detected 11/11 of these samples, demonstrating that the multiplex also detects FMDV in clinical samples that were most likely rendered non-viable for VI between collection and laboratory receipt.

Limit of detection. The limit of detection (LOD) of the 3D and 5'UTR assays in mRT-PCR and rRT-PCR formats were compared using serially diluted clinical samples of serotype O or SAT 2 (Figure 3). At the defined cut-off values, the mRT-PCR LOD was higher than rRT-PCR by approximately 5-625 times for the 3D assay, and 25-125 times for 5'UTR, depending on serotype. The higher LOD of the mRT-PCR therefore caused the loss of diagnostic performance evident in the ROC plots. Further assay optimization efforts are underway to improve the limits of detection of the 3D and 5'UTR assays in the mRT-PCR format.

Differential detection. Characterization of the diagnostic accuracy of the look-alike disease assays within the mRT-PCR assay is ongoing and will be reported separately. The mRT-PCR assay was positive with three SVDV field samples (ITL 4/77, HKN 1/80, HKN 5/91), two VESV isolates

(serotypes -B51 and -H54), San Miguel sea lion virus (serotypes SMSV-7, -9, -10, -11 and -13) and cetacean calicivirus (CCV, Tur-1, dolphin). The mRT-PCR did not detect VESV -B1-34, other VESV serotypes including bovine (Bos-1 (Tillamook)), feline (A4), reptile (rattlesnake), and skunk, and VSV (serotype NJ 15/88 CP211634 and Indiana 1 subtype Ind 2 Maipu Argentina). Feline calicivirus is in a distinct genomic group from the other VESV viruses tested and was not expected to be detected by this assay. mRT-PCR identified look-alike disease viruses in four suspect FMD field samples that had previously been designated FMDV-negative by VI and rRT-PCR (Figure 4). Two samples from cattle (IRN 4/2002, IRQ 58/2002) were mRT-PCR positive for parapoxvirus. The parapox-3 assay generated the strongest response, however all three assays exceeded their respective cut-off values for both samples. Two further cattle samples (UKG 36/94, UKG 37/94) were mRT-PCR positive for BVDV that were also confirmed FMDV-negative by all methods. The mRT-PCR simultaneously detected the co-infection of FMDV and BVDV in a FMDV true positive sample (HUN 2/72). The presence of BVDV, which causes a prevalent disease of cattle, did not mask the detection of FMDV.

Cross-talk between detection channels was minimal, even at the relatively high concentrations of FMDV RNA in many of the field samples tested. A matrix of correlation coefficients for all channel pairings was calculated using the results of all FMDV true positive samples (n=213). The correlation coefficients of the 3D and 5'UTR channels with other channels did not exceed 0.209 and 0.114, respectively.

299

300 **DISCUSSION**

301

The mRT-PCR assay was able to detect FMDV at clinically relevant concentrations. The apparent concentration of FMDV in the clinical samples was relatively high as might be expected with vesicular epithelial tissue which, when infected, is rich in virus. The performance of mRT-PCR and singleplex rRT-PCRs was compared using VI and Ag-ELISA to define FMDV true positives and true negatives.

306 Although VI and Ag-ELISA are established methods for the detection of FMDV, previous studies (14,
307 35) have demonstrated that rRT-PCR has higher diagnostic sensitivity and can detect virus in additional
308 samples, which for the purposes of this study would be classified as “true” negatives. For the majority
309 of true positive field samples, the 3D mRT-PCR response was saturated and grouped far from the cut-
310 off. The 5'UTR mRT-PCR signal was generally lower and clustered on either side of the cut-off. The
311 higher LOD of both assays in the mRT-PCR format may be caused by low-level non-specific
312 interactions between primer sets that could reduce amplification efficiency. The primers used in the
313 5'UTR assay were not originally designed with multiplexing in mind. In order to recognize a wide range
314 of FMDV isolates, the 5'UTR primers have a high degree of degeneracy (32-fold and 8-fold for the
315 forward and reverse primers respectively) that could increase the likelihood of non-specific interactions
316 with other primer sets in the multiplexed reaction mixture. Whilst primer probe sets incorporating
317 degenerate nucleotides can offer broader coverage of highly-variable gene segments, further work is
318 required to refine their design for multiplexed assays. Asymmetric PCR and multivariate optimization
319 may lead to further improvements in the LODs of the 3D and 5'UTR mRT-PCR assays.

320 The 3D and 5'UTR assay responses had serotypic bias, a finding which agrees with earlier
321 observations (4, 14). In mRT-PCR format, the 3D assay was less sensitive for A and C serotypes
322 whereas the 5'UTR was less effective against SAT serotypes. For all serotypes collectively, the 3D
323 assay was more sensitive than the 5'UTR within each format (rRT-PCR 3D 96.7% vs. 5'UTR 87.3%,
324 mRT-PCR 3D 90.1% vs. 5'UTR 58.7%). In an earlier evaluation of the rRT-PCR assays (14), the
325 diagnostic sensitivity of the 3D (97.7 %) was found to be slightly higher than the 5'UTR assay (95.4%).
326 The diagnostic sensitivity increased when the results of the 3D and 5'UTR assays were used in
327 combination. This is due to the “or” nature of the combination, whereby a single mRT-PCR positive
328 result on the 3D or 5'UTR channel generates a combined mRT-PCR positive result. The diagnostic
329 sensitivity of the combined mRT-PCR assay was 93.9%, compared to 98.1% for combined rRT-PCR.
330 The loss of diagnostic sensitivity from rRT-PCR to mRT-PCR, due to higher limits of detection, was

331 partially offset by the inherent ability of the multiplex assay to screen multiple loci simultaneously.
332 Combined rRT-PCR only increased the diagnostic sensitivity by 1.4% compared to using the 3D rRT-
333 PCR assay in isolation which could factor into the cost-benefit of conducting parallel assays. As the
334 mRT-PCR is a screening assay that would most likely be used in conjunction with confirmatory tests,
335 potential users may be more tolerant of lower specificity in order to achieve higher sensitivity. The
336 negative predictive value was also higher for combined assays, because of the “or” nature of the
337 combined result as false negatives only occurred when both the 3D and 5'UTR results agreed.

338 The differential detection of FMDV from look-alike disease viruses, which included ssRNA and
339 dsDNA targets, was demonstrated by testing representative isolates of SVDV, VESV and VSV. For a
340 diagnostic laboratory, this could produce time and cost savings, when compared to testing for each
341 disease using singleplex rRT-PCR assays. For veterinarians, the mRT-PCR assay could increase
342 confidence in a sample identified as FMDV-negative by simultaneously screening for the presence of
343 look-alike diseases. For networks of veterinary diagnostic laboratories, an mRT-PCR assay could
344 facilitate embedded foreign animal disease surveillance whilst conducting routine testing of endemic
345 animal disease viruses. An inter-laboratory evaluation of this multiplex assay was recently conducted in
346 fourteen US National Animal Health Laboratory Network laboratories; the results suggested the mRT-
347 PCR technology could be operated successfully in this setting. The diagnostic performance evaluation
348 for the look-alike disease assays in the mRT-PCR format is underway and will be reported separately.

349 The mRT-PCR format is compatible with procedures and instrumentation used for rRT-PCR. The use
350 of a single method to prepare clinical samples for mRT-PCR and rRT-PCR analysis was demonstrated
351 in this study. The mRT-PCR requires the post-processing of RT-PCR product to the microsphere array
352 which takes ~50 min per 96-well plate using manually operated multi-channel pipettes. The 96-well
353 plate format provides convenient interchangeability between manual and automated platforms. With this
354 reagent set, the Bio-Plex flow cytometer analyzed each well in ~40 s, or ~1 h per 96-well plate.
355 Although the diagnostic sensitivity of the mRT-PCR for FMDV detection is slightly lower than

356 singleplex rRT-PCR, it provides significantly more diagnostic information. With 17 detection channels
357 for seven different viruses, the current prototype panel generates 1632 individual assay results per 96-
358 well plate. The microsphere suspension array is a versatile platform compatible with many different
359 types of diagnostic tests, including immunological and serological assays which could increase its utility
360 within a veterinary diagnostic laboratory, not only for outbreak response and recovery but also for
361 routine testing for endemic diseases. The inherent flexibility of the Luminex array also enables the
362 composition of a given multiplex assay to be altered by simply adding or removing detection channels.
363 Refinements to the first version of the mRT-PCR assay described herein are currently underway,
364 including the development of two species-specific panels for more comprehensive coverage. These new
365 panels incorporate additional assays for other FMDV look-alike diseases. The bovine-specific panel
366 incorporates assays for FMD, Malignant Catarrhal Fever (MCF), Rinderpest, Bluetongue, BHV-1,
367 BVD, Parapox and VS. The porcine-specific panel includes assays for FMD, SVD, VES, VS and
368 Porcine Reproductive and Respiratory Syndrome (PRRS).

369 **ACKNOWLEDGMENTS**

370 This work was performed under the auspices of the U.S. Department of Energy by the University of
371 California Lawrence Livermore National Laboratory under Contract No. W-7405-Eng-48. It was
372 funded by the Department of Homeland Security (DHS) Science and Technology Directorate, award
373 (HSHQDC-06-X-00277) and UK DEFRA project SE1121. The material presented represents the
374 position of the authors and not necessarily that of DHS.

- 376 1. **Baxi, M., D. McRae, S. Baxi, I. Greiser-Wilke, S. Vilcek, K. Amoako, and D. Deregt.** 2006.
377 A one-step multiplex real-time RT-PCR for detection and typing of bovine viral diarrhea viruses.
378 Vet. Microbiol. **116**:37-44.
- 379 2. **Bhudevi, B., and D. Weinstock.** 2001. Fluorogenic RT-PCR assay (TaqMan) for detection and
380 classification of bovine viral diarrhea virus. Vet. Microbiol. **83**:1-10.
- 381 3. **Boyle, D. B., T. Taylor, and M. Cardoso.** 2004. Implementation in Australia of molecular
382 diagnostic techniques for the rapid detection of foot and mouth disease virus. Aust. Vet. J.
383 **82**:421-425.
- 384 4. **Callahan, J. D., F. Brown, F. A. Csorio, J. H. Sur, E. Kramer, G. W. Long, J. Lubroth, S.**
385 **J. Ellis, K. S. Shoulars, K. L. Gaffney, D. L. Rock, and W. M. Nelson.** 2002. Use of a
386 portable real-time reverse transcriptase-polymerase chain reaction assay for rapid detection of
387 foot-and-mouth disease virus. J. Am. Vet. Med. Assoc. **220**:1636-1642.
- 388 5. **Deregt, D., S. A. Gilbert, S. Dudas, J. Pasick, S. Baxi, K. M. Burton, and M. K. Baxi.** 2006.
389 A multiplex DNA suspension microarray for simultaneous detection and differentiation of
390 classical swine fever virus and other pestiviruses. J. Virol. Methods **136**:17-23.
- 391 6. **Dunbar, S. A.** 2006. Applications of Luminex[®] xMAP[™] technology for rapid, high-throughput
392 multiplexed nucleic acid detection. Clin. Chim. Acta **363**:71-82.
- 393 7. **El-Kholy, A. A., S. R. Bolin, J. F. Ridpath, R. M. H. Arab, A. A. Abou-Zeid, H. M.**
394 **Hammam, and K. B. Platt.** 1998. Use of polymerase chain reaction to simultaneously detect
395 and type bovine viral diarrhoea viruses isolated from clinical specimens. Rev. sci. Tech. Off. int.
396 Epiz. **17**:733-742.
- 397 8. **Ferris, N. P., and M. Dawson.** 1988. Routine application of enzyme-linked immunosorbent
398 assay in comparison with complement fixation for the diagnosis of foot-and-mouth and swine
399 vesicular diseases. Vet. Microbiol. **16**:201-209.
- 400 9. **Grubman, M. J., and B. Baxt.** 2004. Foot-and-mouth disease. Clin. Microbiol. Rev. **17**:465-
401 493.
- 402 10. **Han, J., D. C. Swan, S. J. Smith, S. H. Lum, S. E. Sefers, E. R. Unger, and Y.-W. Tang.**
403 2006. Simultaneous amplification and identification of 25 human papillomavirus types with
404 Templex technology. J. Clin. Microbiol. **44**:4157-4162.
- 405 11. **Hole, K., A. Clavijo, and L. A. Pineda.** 2006. Detection and serotype-specific differentiation of
406 vesicular stomatitis virus using a multiplex, real-time, reverse transcription-polymerase chain
407 reaction assay. J. Vet. Diagn. Invest. **18**:139-146.
- 408 12. **James, A. D., and J. Rushton.** 2002. The economics of foot and mouth disease. Rev. sci. Tech.
409 Off. int. Epiz. **21**:637-644.
- 410 13. **Jiménez-Clavero, M. A., M. Agüero, E. S. Miguel, T. Mayoral, M. C. López, M. J. Ruano,**
411 **E. Romero, F. Monaco, A. Polci, G. Savini, and C. Gómez-Tejedor.** 2006. High throughput
412 detection of bluetongue virus by a new real-time fluorogenic reverse transcription-polymerase
413 chain reaction: Application on clinical samples from current Mediterranean outbreaks. J. Vet.
414 Diagn. Invest. **18**:7-17.
- 415 14. **King, D. P., N. P. Ferris, A. E. Shaw, S. M. Reid, G. H. Hutchings, A. C. Giuffre, J. M.**
416 **Robida, J. D. Callahan, W. M. Nelson, and T. R. Beckham.** 2006. Detection of foot-and-
417 mouth disease virus: comparative diagnostic sensitivity of two independent real-time reverse
418 transcription-polymerase chain reaction assays. J. Vet. Diagn. Invest. **18**:93-97.
- 419 15. **Kitching, P. R., P. V. Barnett, D. Paton, D. K. J. Mackay, and A. I. Donaldson.** May 2006,
420 posting date. Foot-and-mouth disease, Manual of Diagnostic Tests and Vaccines for Terrestrial
421 Animals, 5th ed, vol. Part 2, Section 2.1, Chapter 2.1.1. Office International des Epizooties,
422 Paris, France. http://www.oie.int/eng/normes/mmanual/A_00024.htm.
- 423 16. **Kitching, R. P., and G. J. Hughes.** 2002. Clinical variation in foot and mouth disease: sheep

- and goats. *Rev. sci. Tech. Off. int. Epiz.* **21**:505-512.
17. **Knowles, N. J., and A. R. Samuel.** 2003. Molecular epidemiology of foot-and-mouth disease virus. *Virus Res.* **91**:65-80.
 18. **Lee, W.-M., K. Grindle, T. Pappas, D. Marshall, M. Moser, E. Beaty, P. A. Shult, J. Prudent, and J. E. Gern.** 2007. A high-throughput, sensitive and accurate multiplex PCR-microsphere flow cytometry system for large-scale comprehensive detection of respiratory viruses. *J. Clin. Microbiol.* doi:10.1128/JCM.02501-06:JCM.02501-06.
 19. **Létant, S. E., J. I. Ortiz, L. F. Bentley Tammero, J. M. Birch, R. W. Derlet, S. Cohen, D. Manning, and M. T. McBride.** 2007. A multiplexed reverse transcriptase PCR assay for identification of viral respiratory pathogens at point-of-care. *J. Clin. Microbiol.* doi:10.1128/JCM.00880-07:JCM.00880-07.
 20. **Letellier, C., and P. Kerkhofs.** 2003. Real-time PCR for simultaneous detection and genotyping of bovine viral diarrhea virus. *J. Virol. Methods* **114**:21-27.
 21. **Li, H., M. A. McCormac, R. W. Estes, S. E. Sefers, R. K. Dare, J. D. Chappell, D. D. Erdman, P. F. Wright, and Y.-W. Tang.** 2007. Simultaneous detection and high-throughput identification of a panel of RNA viruses causing respiratory tract infections. *J. Clin. Microbiol.* **45**:2105-2109.
 22. **Mahlum, C. E., S. Haugerud, J. L. Shivers, K. D. Rossow, S. M. Goyal, J. E. Collins, and K. S. Faaberg.** 2002. Detection of bovine viral diarrhea virus by TaqMan[®] reverse transcription polymerase chain reaction. *J. Vet. Diagn. Invest.* **14**:120-125.
 23. **Mahony, J., S. Chong, F. Merante, S. Yaghoubian, T. Sinha, C. Lisle, and R. Janeczko.** 2007. Development of a Respiratory Virus Panel (RVP) Test for the Detection of Twenty Human Respiratory Viruses using Multiplex PCR and a Fluid Microbead-based Assay. *J. Clin. Microbiol.* doi:10.1128/JCM.02436-06.
 24. **NCCLS.** 1995. Assessment of the clinical accuracy of laboratory tests using receiver operating characteristic (ROC) plots. NCCLS, Wayne, PA.
 25. **NCCLS.** 2002. User protocol for evaluation of qualitative test performance; approved guideline. EP12-A. NCCLS, Wayne, PA.
 26. **Nitsche, A., M. Buttner, S. Wilhelm, G. Pauli, and H. Meyer.** 2006. Real-time PCR detection of parapoxvirus DNA. *Clin. Chem.* **52**:316-319.
 27. **Orrú, G., M. L. Ferrando, M. Meloni, M. Liciardi, G. Savini, and P. De Santis.** 2006. Rapid detection and quantitation of Bluetongue virus (BTV) using a Molecular Beacon fluorescent probe assay. *J. Virol. Methods* **137**:34-42.
 28. **Rasmussen, T. B., A. Uttenthal, and M. Agüero.** 2006. Detection of three porcine vesicular viruses using multiplex real-time primer-probe energy transfer. *J. Virol. Methods* **134**:176-182.
 29. **Reid, S. M., N. P. Ferris, G. H. Hutchings, D. P. King, and S. Alexandersen.** 2004. Evaluation of real-time reverse transcription polymerase chain reaction assays for the detection of swine vesicular disease virus. *J. Virol. Methods* **116**:169-176.
 30. **Reid, S. M., N. P. Ferris, G. H. Hutchings, Z. D. Zhang, G. J. Belsham, and S. Alexandersen.** 2002. Detection of all seven serotypes of foot-and-mouth disease virus by real-time, fluorogenic reverse transcription polymerase chain reaction assay. *J. Virol. Methods* **105**:67-80.
 31. **Reid, S. M., D. P. King, A. E. Shaw, N. J. Knowles, G. H. Hutchings, E. J. Cooper, A. W. Smith, and N. P. Ferris.** 2007. Development of a real-time reverse transcription polymerase chain reaction assay for detection of marine caliciviruses (genus Vesivirus). *J. Virol. Methods* **140**:166-173.
 32. **Schmitt, M., I. G. Bravo, P. J. F. Snijders, L. Gissmann, M. Pawlita, and T. Waterboer.** 2006. Bead-based multiplex genotyping of human papillomaviruses. *J. Clin. Microbiol.* **44**:504-512.
 33. **Shaw, A. E., P. Monaghan, H. O. Alpar, S. Anthony, K. E. Darpel, C. A. Batten, A. Guercio, G. Alimena, M. Vitale, K. Bankowska, S. Carpenter, H. Jones, C. A. L. Oura, D.**

- P. King, H. Elliott, P. S. Mellor, and P. P. C. Mertens.** 2007. Development and initial evaluation of a real-time RT-PCR assay to detect bluetongue virus genome segment 1. *J. Virol. Methods* doi:10.1016/j.jviromet.2007.05.014.
34. **Shaw, A. E., S. M. Reid, K. Ebert, G. H. Hutchings, N. P. Ferris, and D. P. King.** 2007. Implementation of a one-step real-time RT-PCR protocol for diagnosis of foot-and-mouth disease. *J. Virol. Methods* **143**:81-85.
35. **Shaw, A. E., S. M. Reid, D. P. King, G. H. Hutchings, and N. R. Ferris.** 2004. Enhanced laboratory diagnosis of foot and mouth disease by real-time polymerase chain reaction. *Rev. sci. Tech. Off. int. Epiz.* **23**:1003-1009.
36. **Slezak, T., T. Kuczmarski, L. Ott, C. Torres, D. Medeiros, J. Smith, B. Truitt, N. Mulakken, M. Lam, E. Vitalis, A. Zemla, C. Zhou, and S. Gardner.** 2003. Comparative genomics tools applied to bioterrorism defence. *Brief Bioinform* **4**:133-49.
37. **Traul, D. L., S. Elias, N. S. Taus, L. M. Herrmann, J. L. Oaks, and H. Li.** 2005. A real-time PCR assay for measuring alcelaphine herpesvirus-1 DNA. *J. Virol. Methods* **129**:186-190.
38. **Vignali, D. A. A.** 2000. Multiplexed particle-based flow cytometric assays. *J. Immunol. Methods* **243**:243-255.
39. **Young, N. J., C. J. Thomas, M. E. Collins, and J. Brownlie.** 2006. Real-time RT-PCR detection of bovine viral diarrhoea virus in whole blood using an external RNA reference. *J. Virol. Methods* **138**:218-222.
40. **Zou, S., J. Han, L. Wen, Y. Liu, K. Cronin, S. H. Lum, L. Gao, J. Dong, Y. Zhang, Y. Guo, and Y. Shu.** 2007. Human influenza A virus (H5N1) detection by a novel multiplex PCR typing method. *J. Clin. Microbiol.* **45**:1889-1892.

498 **FIGURE CAPTIONS**

499

500 **FIG. 1.** Schematic of the multiplex RT-PCR assay. In the presence of target nucleic acid, the
501 biotinylated forward primer is extended during the PCR. The PCR product is hybridized to the
502 microsphere array, whereby the extended forward primer binds to the complementary probe-labeled
503 microsphere. The complex is labeled with fluorescent reporter (SAPE) then analyzed using a Bio-Plex
504 flow cytometer. The fluorescence (570 nm) of bound reporter molecules is measured, and the median
505 fluorescence intensity is calculated for each microsphere class, then compared to a cut-off to indicate the
506 presence or absence of a target nucleic acid sequence in a sample.

507

508 **FIG. 2.** Receiver Operating Characteristic (ROC) plots for the RNA polymerase (3D) and the
509 5'untranslated region (5'UTR) assays in singleplex real-time reverse transcriptase PCR (rRT-PCR) and
510 multiplexed reverse transcriptase polymerase chain reaction (mRT-PCR) formats. True positive
511 samples ($n=213$) representative of all FMDV serotypes and true negative samples ($n=74$) were
512 analyzed. The area under the curves are; 3D rRT-PCR (0.985), 5'UTR rRT-PCR (0.942), 3D mRT-
513 PCR (0.955) and 5'UTR mRT-PCR (0.773).

514

515 **FIG. 3.** Comparison of the RNA polymerase (3D) and 5'untranslated region (5'UTR) assay responses
516 in multiplex RT-PCR and singleplex rRT-PCR formats by consecutive 5-fold serial dilutions of three
517 different FMDV true positive epithelial suspension samples (undiluted to more dilute from right to left).
518 The multiplex RT-PCR detection was less sensitive than singleplex rRT-PCR (approximately 5-625
519 times and 25-125 times, respectively for 3D and 5'UTR signatures). Cut-off values for rRT-PCR ($C_t \leq$
520 32) and mRT-PCR ($3D\ MFI \geq 6.5$, $5'UTR \geq 5.5$) are indicated by the vertical and horizontal dashed
521 lines, respectively. The 5'UTR signature did not respond to the SAT 2 (SAU 4/2000) sample at any

522 dilution in singleplex or multiplex formats and was omitted from the plot. Error bars indicate $\pm 1\sigma$ of the
523 mean (n=2) response from the signature in mRT-PCR format. For rRT-PCR each sample was analyzed
524 in singlet (n=1).

525

526 **FIG. 4.** Multiplex RT-PCR assay identification of FMDV look-alike disease viruses in suspect FMDV
527 clinical sample submissions. The mRT-PCR assay ruled out FMDV whilst simultaneously ruling in
528 FMD look-alike disease viruses. A and B show the multi-loci detection of parapox viruses. C shows
529 BVDV detection in field samples from cattle (UKG 36/94, UKG 37/94) that were tested FMDV-
530 negative by virus isolation and rRT-PCR and the simultaneous detection of FMDV and BVDV (HUN
531 2/72). mRT-PCR cut-offs indicated by the vertical dashed lines were PPOX 1 (≥ 7.5), PPOX 3 (≥ 9.5)
532 and BVDV (≥ 6.5). The horizontal dashed lines indicate the cut-off for FMDV 3D (≥ 6.5).

533 **TABLE 1.** Primer and probe sequences of the multiplex RT-PCR assay.

Assay name	Forward Primer (5'→3') Reverse Primer (5'→3')	Probe (5'→3')
BHV 1	GTGCCAGCCGCGFAAAAG GACGACTCCGGGCTCTTTT	TCCTGGTTCCAGAGCGCTAACATGGAG
BHV 2	TGAGGCCFATGTATGGGCAGFT GCGCGCCAAACATAAGTAAA	AAATAACACGGTGTGCACTTAAATAAGATTTCGG
BTB 1	GCACCCFATATGTTTCCAGACCA CAGCTAACTCTTCAGCCACACG	CTAACTCGTGGGCCAATCATCATCTTCTGT
BTB 2	AGAATTCAGGAFTGGGCAGGA GCACAATTCCTATCCCCTTA	CCATCACACCATTATACTGTACCCGCGTAGC
BVDV	GGTAGTCGTCAGTGGTTCGAC CATGTGCCATGTACAGCAGAGAT	CCTCGTCCACGTGGCATCTCGAG
FMDV 3D	ACTGGGFTTTACAAACCTGTGA GCGAGTCCTGCCACGGA	GTCCCACGGCGTGCAAAGGA
FMDV 5'UTR	CACYTYAAGRFGACAYTGRTACTGGTAC CAGATYCCRAGTGWCICITGTTA	CCTCGGGGTACCTGAAGGGCATCC
PPOX 1	GCAGATGCGCTCCTGGTT GCACCTCTGCTGCTGCAA	CCGACTCCGACGTGGAGAACGTG
PPOX 2	GATGGCCGTCAGCTCTT CGTACAAGATCACGGCCAACT	TGTACGGGCTCATGGGCTTCCG
PPOX 3	GCAGCAGTGCACCACGFACT CGCTGAACCCGTACATCCT	GACTTCGAGGCGGACAACAAGCG
SVDV 1	CAGGAFAATTTCTFCCAAGGGC ACGTGAACATTTTCGAGCTTCC	TGCATTGTGTCTGATGGTACAACCTGTGACG
SVDV 2	GACTTGFTGTGGCTGGAGGA CAGCGCCATGGTGAGGTAG	TGACCGTAATGAGGTCATCGTGATTTCTCAC
SVDV 3	GACAAAGTGGCCAAGGGAAA CACGTAAACCACACTGGGCT	CTGGCGTCATAGCCTGAATAGTCAAACGCTA
VESV 1	GCCTFCTCCCTFCCCAAAA TGAAGGAATGGTTCCGTCAGT	CATCATCGTTGATAACCTTAGATGTGCAATTTGG
VESV 2	GGGAATFAGGGTGTGCAFCATT CACGTCTTGATGTTGGCTTGAC	AAATTGGCATAATCAACCTTGTCAGATGAGTCG
VESV 3	GGTCGCTCTCACTGATGATGAGTA GGTGTATATCAGCACCCATTGC	GCTCGGTGCCTGAGTTGGAGGAAG
VESV 4	ACCACCTCTGGAAACATCFATGG TTTGTGCACGTGTCACGAAT	CGGGACGGGCATTTGTCACCA
FC	N/A	CAAAGTGGGAGACGTCGFTG
IC	N/A	CAAAGTGGGAGACGTCGTTG-Cy3
NC	N/A	CAAAGTGGGAGACGTCGTTG

534

535 BHV = bovine herpes virus-1, PPOX = parapox virus complex, FMDV = foot-and-mouth disease virus,
536 BVDV = Bovine Viral Diarrhea Virus, BTB = bluetongue virus, SVDV = swine vesicular disease virus,
537 VESV = vesicular exanthema of swine virus, NC = negative control, FC = fluorescence control, IC =
538 instrument control, Cy3 = fluorescent cyanine dye, Y = pyrimidines (C/T), R = purines (A/G), W =
539 weak 2-bonds (A/T), I = inosine (universal base). F = an internal biotinylated dT. All forward primers
540 also include biotinylation at the 5' terminus. All probes contain an amine attached to the 5' terminus
541 with a carbon 6 and internal spacer 18. N/A = not applicable.

542

543 **TABLE 2.** Performance metrics for the 3D and 5'UTR FMDV assays in rRT-PCR and mRT-PCR
544 formats using independent or combined results.

Performance metric		rRT-PCR			mRT-PCR		
		3D	5'UTR	Combined	3D	5'UTR	Combined
Sensitivity (%) (Fraction)	A	93.5 (43/46)	87.0 (40/46)	97.8 (45/46)	80.4 (37/46)	76.1 (35/46)	93.5 (43/46)
	Asia 1	100 (10/10)	100 (10/10)	100 (10/10)	100.0 (10/10)	100 (10/10)	100 (10/10)
	C	95.0 (19/20)	95.0 (19/20)	95.0 (19/20)	85.0 (17/20)	55.0 (11/20)	90.0 (18/20)
	O	96.7 (87/90)	86.7 (78/90)	97.8 (88/90)	91.1 (82/90)	70.0 (63/90)	92.2 (83/90)
	SAT 1	100.0 (18/18)	88.9 (16/18)	100 (18/18)	100 (18/18)	16.7 (3/18)	100 (18/18)
	SAT 2	100.0 (21/21)	71.4 (15/21)	100 (21/21)	95.2 (20/21)	9.5 (2/21)	95.2 (20/21)
	SAT 3	100.0 (8/8)	100 (8/8)	100 (8/8)	100 (8/8)	12.5 (1/8)	100 (8/8)
	All	96.7 (206/213)	87.3 (186/213)	98.1 (209/213)	90.1 (192/213)	58.7 (125/213)	93.9 (200/213)
	C.I. (95%)	93.4-98.4	82.2-91.1	95.3-99.3	85.4-93.5	52-65.1	89.8-96.4
	Specificity (%)	94.6	95.9	93.2	93.2	94.6	91.9
C.I. (95%)		86.9-97.9	88.7-98.6	85.1-97.1	85.1-97.1	86.9-97.9	83.4-96.2
Sample Disease Prevalence (%)		74.2	74.2	74.2	74.2	74.2	74.2
PPV (%)		98.1	98.4	97.7	97.5	96.9	97.1
NPV(%)		90.9	72.4	94.5	76.7	44.3	84.0
Efficiency (%)		96.2	89.5	96.9	90.9	67.9	93.4

545 Sensitivity is the fraction of the 213 true positive epithelial tissue samples (subjected to both multiplex

546 RT-PCR and rRT-PCR assays) that are assigned as FMDV-positive by the indicated assay. C.I.
 547 Confidence Interval, PPV Positive Predictive Value, NPV Negative Predictive Value.
 548

549 **TABLE 3.** 3-Way comparison table for the combined results of 3D and 5'UTR assays in rRT-PCR and
 550 mRT-PCR formats.

Method Result		Total Samples	True Diagnosis	
Combined	Combined		Positive	Negative
rRT-PCR	mRT-PCR			
Positive	Positive	202	200	2
Positive	Negative	12	9	3
Negative	Positive	4	0	4
Negative	Negative	69	4	65
Total		287	213	74

551

FIG.1




FIG.2




FIG. 3





FIG. 4

